

Imaging and Resection

Next Generation in Minimally Invasive Surgical Technology



Arthrex® 

Helping Surgeons Treat Their Patients Better[®]

Since its inception, Arthrex has been committed to one mission: Helping Surgeons Treat Their Patients Better. We are strategically focused on constant product innovation through scientific research, surgeon collaboration, and medical education to make less invasive surgical procedures simple, safer, and more reproducible. Each year, we develop more than 1,000 new innovative products and procedures to advance minimally invasive orthopedics worldwide.

Arthrex has always remained a privately held company, which allows for the rapid evaluation of new technologies and ideas and the freedom to develop products and techniques that truly make a difference. Our experienced team of dedicated professionals represents a shared passion and commitment to delivering uncompromising quality to the health care providers who use our products and the millions of patients whose lives we impact.

The medical significance of our contributions serves as our primary benchmark of success and will continue into the future as the legacy of Arthrex.


05	Integration, Imaging, and Image Management
27	Nano Arthroscopy
45	Synergy Power™ System
57	Synergy ^{Resection} ™ System
67	Fluid and Insufflation Management
81	Arthroscopic RF System
89	Synergy Service and Support

Integration, Imaging, and Image Management

06	Synergy Matrix™ Video Integration
07	Synergy.net™ Data Integration Software
07	SurgeonVault® Application
08	Synergy OR Command™ Room Status
09	Synergy 4K Imaging Platforms
11	UHD 4K Monitors
11	4K and HD Camera Heads
12	Light Guides
12	StarFish™ Scope Rotation Attachment
14	Arthroscopes and Sheath Systems
17	Laparoscopes
17	HD Endoscopes
18	Scope and Sheath Accessories
18	Competitive Sheaths
18	Instrument Cases
19	Bridges
19	Metal Cannulas and Sets
21	Video Cart
22	Video Cart and Monitor Stand
24	Secondary Monitor Stand
24	Printer and Recorder

Integration

Synergy Matrix™ Video Integration



The Synergy Matrix system was the world's first OR integration platform capable of supporting the 4K video signal produced by the Synergy Vision™ and Synergy^{ID}™ camera systems and distributing it to multiple displays within the OR. The next-generation platform expands on the success of the Synergy Matrix integration system with the inclusion of an integrated touch panel that provides improved system performance, reliability, and overall user experience.

Please direct questions or requests for customized Synergy Matrix integration quotes to the Synergy Integration team at SynergyIntegration@Arthrex.com.

Synergy Matrix software license	AR-3200-1052
Synergy Matrix Core software license	AR-3200-1053

Key Features and Benefits

- › A redesigned user interface is available on a touchscreen in every integrated OR.
- › Programmable room presets allow for quick routing and configuration of a room's audio and video sources to destinations. As an example, a room preset can be set for laparoscopic procedures so with the touch of a button, the room automatically routes the configured sources to the desired displays in the room.
- › Robust room-to-room communication and source sharing gives users more control over which outside rooms have viewing capabilities of their internal sources. Shared sources can be granted permission for a specific time range or removed from sharing by the source owner as needed.

- › Room camera control allows users to modify the pan, tilt, and zoom settings for the room camera. In addition, room camera settings can be saved as a preset to allow for quick configuration of a room camera's position. For example, a room camera preset can be set for privacy so with the touch of a button, the room camera's position automatically adjusts away from the patient's location.
- › Improved remote support and proactive system monitoring allow for integrated rooms to be upgraded and serviced independently.
- › Equipped with the latest in 4K imaging and integration technology, the Synergy Matrix system is a state-of-the-art, video-over-IP solution that can lower the cost of ownership and positively impact OR workflow.
- › Room-to-room connectivity supports sharing video and audio sources across integrated ORs, conference rooms, pathology labs, and observation rooms.
- › The Synergy Matrix system provides more capabilities at a fraction of the footprint required in the OR, thereby saving valuable OR space.

Synergy Matrix Core

The Synergy Matrix Core system provides in-room 4K video and audio integration while simplifying the installation process and focusing on key functionality. The value-driven platform uses the same touch panel controller as the Synergy Matrix system but with a simplified user interface, providing both uniformity and familiarity for the OR staff, while featuring intuitive and easy-to-use functionality.

Additionally, an all-in-one switch assembly installed in each integrated room allows for true 4K video and audio routing, operating completely independently of the facility's network and without cable pulls to a centralized IT closet.

Synergy.net™ Data Integration Software



Synergy.net data integration software expands the network capabilities of the Synergy Vision™ surgical video system by synchronizing surgeon and procedure settings across each connected OR to provide a consistent and reliable user experience. In addition, improved OR workflows can be achieved by integrating with your electronic health record (EHR), PACS, surgical schedule, and vendor neutral archive.

A lightweight software-based application, the Synergy.net system is installed on a facility-provided virtual machine. This approach gives the facility's IT full control of data activities, security, and HIPAA requirements. Automated system health monitoring services are leveraged to provide proactive system support.

For surgeons with premium SurgeonVault accounts, the Synergy.net system can be configured to upload case media automatically to the SurgeonVault secure cloud-based surgical communication application.

Synergy.net HL7 bidirectional interface	AR-3200-1034
Synergy.net HL7 single interface	AR-3200-1047
Synergy.net HL7 additional facility	AR-3200-1055
Synergy.net pairing license key	AR-3200-1067
Synergy.net network printing	AR-3200-1043

- › Interfaces with EHRs and other HL7-based systems
- › Provides a single database for all cases captured by Synergy devices
- › Improves DICOM capabilities, which includes transfer of video to a PACS
- › Offers network printing from any Synergy device

SurgeonVault® Application




The SurgeonVault secure cloud-based surgical communication application provides surgeons with access to their data at any time and allows them to distribute surgical videos and stills and other selected content directly to their patients.

Key Features and Benefits

- › Improved surgeon-patient communication
 - › Surgeons can quickly and easily share case and library content to meet their patients' needs
 - › Patients receive access to a secure, multimedia-rich patient portal rather than traditional static content
 - › Improved communication to patients may improve patient satisfaction scores
- › User-friendly interface with templates and library content for quick, easy surgeon onboarding
- › Social media tools to promote practice online presence
- › Surgeons can access content from either the Synergy Surgeon™ app or from any web browser at any time from any location
- › Latest encryption technology ensures content remains private and secure for patients and surgeons

SurgeonVault single-user license, 1-year term	SVSIN-1YR
SurgeonVault enterprise license unlimited users, 1-year term	SVENT-1YR

Synergy OR Command™ Room Status



A powerful tool for maximizing OR efficiency and workflow, Synergy OR Command room status provides surgeons and staff with a centralized dashboard view of simultaneous video streams across all operating rooms. Authorized users can access Synergy OR Command room status through a network browser or smart device.

Benefits of Synergy OR Command room status include adaptive scheduling, minimized downtime, accelerated surgeon workflows, and real-time information for decision-makers. Whether viewing the user interface on a computer, tablet, or dedicated wall display in an OR corridor, Synergy OR Command room status provides the efficiency needed in today's busy OR.

Key Features and Benefits

- › Improve workflow efficiencies with OR visibility
- › Kiosk Mode for dedicated displays
- › Privacy Mode allows staff to mute video feeds
- › Configurable dashboard layout options
- › OR turnover status
- › Room camera control
- › Scope view from Synergy camera consoles

Synergy OR Command room status

AR-3200-1044

Imaging and Imaging Management

Synergy 4K Imaging Platforms



Experience the Synergy Vision™ imaging system, an all-in-one option featuring high dynamic range (HDR) 4K visualization with wide-view imaging and fluorescence while supporting the NanoNeedle™ scope and video integration.

Key Features and Benefits

The Synergy Vision imaging system was designed with the latest technology to deliver optimal and revolutionary results.

- › **Pano™ Scope:** Quickly and easily switch from traditional 30° and 70° views to a 45° ultrawide view at the touch of a button. The Pano scope allows surgeons to visualize more anatomy on a single screen while minimizing the need for additional portals or scopes.
- › **Nano Vision™ Functionality:** The NanoNeedle scope integrates directly to the Synergy Vision console as either a primary or secondary view allowing for simultaneous 4K and Nano visualization on the same screen.
- › **Synergy Vision Connect™ Console:** With 6 inputs and 4 outputs, the Synergy Vision Connect system offers built-in OR integration capabilities for in-room switching and routing.
- › **HDR:** The Synergy Vision imaging system includes HDR for enhanced contrast of surgical anatomy and ideal visualization.
- › **Fluorescence Imaging:** Switch to fluorescence 4K imaging with the touch of a button. Easily toggle through different modes and colors, depending on the fluorescence application and individual surgeon preference.

Synergy Tablet Controller

Our unique and innovative tablet interface allows users to remotely and efficiently enter patient information, surgeon preferences, and image management settings once. The Synergy tablet controller may help reduce cost, shorten OR times, and increase efficiency.

Synergy 4K Consoles

Synergy ^{HD} ™ 4K console	AR-3200-0025
Synergy Vision console	AR-3200-0026
Synergy Vision Connect console (a)	AR-3200-0027
Synergy Vision tablet (b)	AR-3200-1016
Synergy digital documentation tablet	AR-3200-1021

The camera console comes with base software to manage 4K image and video capture, printing, and USB export. Additional software modules are offered to better address the customer's individual needs. By offering software modules focused on functionality, facilities can select the features that best serve their unique needs.

Synergy Software Modules

Premium bundle	AR-3200-1091
Data management module	AR-3200-1092
Workflow support module	AR-3200-1093
Clinical enhancement module	AR-3200-1094

Base Module Features

- › Still image capture and video recording
- › Print still images
- › 4K image display
- › Export to USB

Data Management Module Features

- › Export media to DICOM or network folders
- › Synergy.net™ pairing license

Workflow Support Module Features

- › Network-based system allows live video streaming to any authorized remote viewer
- › Export images to the Synergy Surgeon™ app:
 - › Review and edit still images and videos
 - › Annotate on stills
 - › Create presentations and postoperative reports for patients
- › Heads-up display
- › Capture images from a secondary device
- › Surgical preparation checklist
- › 4K video capture

Clinical Enhancement Module Features

- › Augmented reality overlay to enhance spectral and textural changes in the surgical image
- › Enhance or reduce the color red in live image to improve visualization
- › Consistent illumination across the entire image with highlight feature
- › Pano™ scope wide-view imaging

Premium Bundle

- › Comprehensive bundle that includes all Synergy software features, except the Synergy Vision Connect license

Synergy Vision Connect™ System

Synergy Vision Connect software license	AR-3200-1088
---	--------------

Synergy Accessories

Synergy System Integration Cable Kit	AR-3200-1040
Input/output 1080p video signal converter	AR-3210-1005
C-mount HD optical zoom coupler	AR-3210-0006
C-mount right angle coupler	AR-3210-0009
C-mount optical non-zoom coupler	AR-3210-0012
4K wireless receiver	LBW-5-H2-12G-R-A
4K wireless transmitter	LBW-5-H2-12G-T-A
Synergy laser light source	AR-3200-1018

UHD 4K Monitors



- › 27 in, 32 in, and 55 in UHD 4K display options available
- › Versatile multiformat signal support
- › With and without Synergy Matrix™ video integration enabled

Monitors

32 in 4K HDR 12G SDI display	AR-3250-3212
32 in 4K surgical display	AR-3250-3211
27 in 4K surgical display	AR-3250-2701
55 in 4K medical grade display w/ decoder	K9352362
55 in 4K medical grade display 12G version	K9352361
55 in 4K medical grade display non-MNA version	K9352360

Monitor Covers

Protection cover for 32 in monitor, black	703459-15
Protection cover for 27 in monitor, black	703458-15
Tablet cover for Synergy ^{UHD4} ™ system, black	703510-15

4K and HD Camera Heads




The Synergy Vision™ imaging system supports two camera heads. With unique features and benefits for each style, these camera heads were designed for both orthopedic and endoscopic procedures. The picture produced by all of the Synergy camera heads increases depth perception and creates more precise color reproduction. All of these camera heads carry the industry-leading standard autoclavable warranty while maintaining a comfortable and ergonomic design.

Camera Heads

4K Synergy ^{UHD4} camera head, broadband, autoclavable	AR-3210-0029
4K Synergy ^{UHD4} C-mount camera head, broadband, autoclavable	AR-3210-0030
4K Synergy ^{UHD4} ultra camera head, autoclavable	AR-3210-0031
4K Synergy ^{UHD4} ultra C-mount camera head, autoclavable	AR-3210-0032
Synergy ^{HD} 4MOS camera head	AR-3210-0033
Synergy Visionary™ camera head, 4K, C-mount, autoclavable	AR-3210-0034
Synergy Vision camera head, 4K, C-mount, autoclavable	AR-3210-0035
Synergy Vision C-mount camera head, 0°, autoclavable	AR-3210-0036

Light Guides



Light guides are made from hundreds of high-quality fibers. Instead of gluing fibers, we fuse ours at the proximal end to maximize light transmission. The unique handle increases durability and prevents accidental removal of the light guide. Additionally, our transparent casing makes it easier to monitor when to replace light guides.

Light guide cable, blue, 5 mm × 2.7 m (a)	AR-3240-5027
Light guide cable, grey, 5 mm × 2.7 m	AR-3240-5028
Light guide cable, blue, 3.5 mm × 2.7 m	AR-3240-3527
Light guide cable, blue, 5 mm × 4 m	AR-3240-5040
Light guide cable, grey, 5 mm × 2.7 m, 90° distal end	SS9-GYW50W108
Endoscope cable, gray, 3.5 mm	SS-GBS35S120
Light guide cable, grey, universal, 3.5 mm × 2.7 m	AR-3240-3530
Bifurcated light guide cable, grey, 3.5 mm × 2.7 m	AR-3240-5045
Bifurcated light guide cable, gray, 90° distal end, 3.5 mm × 2.7 m	AR-3240-5090
Light guide, blue, 5 mm × 2.7 mm (b)	FTIMG26203TB
Light guide, blue, 3.5 mm × 2.7 mm	FTIMG26206TB

Adapters for Universal Light Guide

Wolf adapter	AR-3240-3530W
Storz adapter	AR-3240-3530S
Olympus adapter	AR-3240-3530O
ACMI adapter	AR-3240-3530A

StarFish™ Scope Rotation Attachment




The StarFish scope rotation attachment connects to Synergy 4 mm and 3 mm arthroscopes and 5.5 mm laparoscopes. The StarFish attachment protrudes proximally toward the camera head, which allows one-handed directional changes without the use of the light cord.

StarFish scope attachment	AR-3370-0006
StarFish C-mount scope attachment	AR-3370-0008



Arthroscopes and Sheath Systems



4K Pano™ Scope

4K C-mount Pano arthroscope, 45°, 4.2 mm × 161 mm (a)	AR-3355-4050
---	--------------

Sheaths and Obturators for Pano Scope

High-flow sheath for Pano scope, tapered, fenestrated, single stopcock	AR-3373-4201
High-flow sheath for Pano scope, tapered, fenestrated, double stopcock (b)	AR-3373-4202
Conical obturator for Pano scope with handle	AR-3375-4201
Sharp trocar for Pano scope with handle	AR-3375-4202
Cannulated obturator for Pano scope with handle	AR-3375-4203
Blunt obturator for Pano scope with handle	AR-3375-4204
Conical obturator for Pano scope	AR-3375-4207
Switching stick, 4.2 mm × 15 in	AR-3020-10
Switching stick, cannulated, 4.2 mm × 15 in	AR-3020-11
Sheath, Hi-flow, 1 stopcock for Pano scope	AR-3372-4201
Sheath, Hi-flow, 2 stopcock for Pano scope	AR-3372-4202
Conical Obturator, for hi-flow Pano scope sheath, with handle	AR-3375-4211
Trocar for hi-flow Pano scope sheath, with handle	AR-3375-4212
Cannulated obturator for hi-flow Pano scope sheath, with handle	AR-3375-4213
Blunt obturator for hi-flow Pano scope sheath, with handle	AR-3375-4214
Conical obturator, for hi-flow pano scope sheath	AR-3375-4217

Inflow/Outflow Sheath for Pano Scope

Inflow/outflow sheath set for Pano scope	AR-3370-4203
Inner sheath for Pano scope, part A	AR-3370-4203A
Outer sheath for Pano scope, part B	AR-3370-4203B
Conical obturator for inflow/outflow sheath for Pano scope	AR-3375-4201A
Sharp trocar for inflow/outflow sheath for Pano scope	AR-3375-4202A

4.8 mm 4K Arthroscopes

4K arthroscope, 30°, 4.8 mm × 152.5 mm	AR-3350-5230
4K arthroscope, 70°, 4.8 mm × 156 mm	AR-3350-5270
4K C-mount arthroscope, 30°, 4.8 mm × 152.5 mm	AR-3355-5230
4K C-mount arthroscope, 70°, 4.8 mm × 156 mm	AR-3355-5270

Sheaths and Obturators for 4.8 mm 4K Arthroscopes

Sheath for 4.8 mm scopes, single stopcock	AR-3373-5001
Sheath for 4.8 mm scopes, double stopcock	AR-3373-5002
Conical obturator w/ handle (c)	AR-3375-5001
Sharp trocar w/ handle	AR-3375-5002
Cannulated obturator w/ handle	AR-3375-5003
Conical obturator (d)	AR-3375-5007
Switching stick, 4.8 mm × 38 cm	AR-3020-05

4 mm 4K Arthroscopes

4K arthroscope, 30°, 4 mm × 152.5 mm	AR-3350-4031
4K arthroscope, 70°, 4 mm × 156 mm	AR-3350-4071
4K C-mount arthroscope, 30°, 4 mm × 152.5 mm	AR-3355-4031
4K C-mount arthroscope, 70°, 4 mm × 156 mm	AR-3355-4071

4 mm Arthroscopes HD

HD arthroscope, 30°, 4 mm × 152.5 mm	AR-3350-4030
HD arthroscope, 30°, 4 mm × 152.5 mm, reverse light post	AR-3350-4030R
C-mount HD arthroscope, 0°, 4 mm × 152.5 mm	AR-3355-4000
C-mount HD arthroscope, 30°, 4 mm × 152.5 mm, reverse light post	AR-3355-4030R

Sheaths and Obturators for 4 mm Arthroscopes

Sheath for 4 mm scopes, single stopcock	AR-3372-4001
Sheath for 4 mm scopes, double stopcock	AR-3372-4002
Conical obturator w/ handle	AR-3375-4011
Sharp trocar w/ handle	AR-3375-4012
Conical obturator	AR-3375-4017
Switching stick, 4.6 mm × 38 cm	AR-3020-07
High-flow sheath for 4 mm scopes, tapered, fenestrated, single stopcock	AR-3373-4001
High-flow sheath for 4 mm scopes, tapered, fenestrated, double stopcock (b)	AR-3373-4002
Conical obturator w/ handle	AR-3375-4001
Sharp trocar w/ handle	AR-3375-4002
Cannulated obturator w/ handle	AR-3375-4003
Blunt obturator w/ handle	AR-3375-4004
Conical obturator	AR-3375-4007
Switching stick, 4 mm × 23 cm	AR-3020
Switching stick, 4 mm × 31 cm	AR-3026
Switching stick, 4 mm × 38 cm	AR-3020-01
Switching stick, 4 mm × 46 cm	AR-3020-02
Switching stick, cannulated, 4 mm × 38 cm	AR-6520-15
Sheath for 4 mm scopes, fenestrated, single stopcock	AR-3374-4001
Sheath for 4 mm scopes, fenestrated, double stopcock	AR-3374-4002
Conical obturator w/ handle	AR-3374-4021
Sharp trocar w/ handle	AR-3375-4022
Conical obturator	AR-3375-4027

Inflow/Outflow Sheath for 4 mm Arthroscopes

Sheath system, in-/outflow 2 stopcock for 4 mm scope	AR-3370-4003
Inner sheath for inflow/outflow Synergy 4 mm scope sheath	AR-3370-4003A
Outer sheath for inflow/outflow Synergy 4 mm scope sheath	AR-3370-4003B
Conical obturator	AR-3375-4005

3 mm Arthroscopes HD

HD arthroscope, 30°, 3 mm × 138 mm	AR-3350-3030
HD arthroscope, 70°, 3 mm × 140 mm	AR-3350-3070
C-mount HD arthroscope, 30°, 3 mm × 138 mm	AR-3355-3030
C-mount HD arthroscope, 70°, 3 mm × 140 mm	AR-3355-3070

Sheaths and Obturators for 3 mm Arthroscopes

Sheath for 3 mm scopes, single stopcock	AR-3372-3001
Sheath for 3 mm scopes, double stopcock	AR-3372-3002
Conical obturator w/ handle	AR-3375-3001
Sharp trocar w/ handle	AR-3375-3002
Conical obturator	AR-3375-3007
High-flow sheath for 3 mm scopes, single stopcock	AR-3372-0301
High-flow sheath for 3 mm scopes, double stopcock	AR-3372-0302
Conical obturator w/ handle	AR-3375-0311
Sharp trocar w/ handle	AR-3375-0312
Conical obturator	AR-3375-0317
High-flow sheath for 3 mm scopes, tapered, fenestrated, single stopcock	AR-3372-0301
High-flow sheath for 3 mm scopes, tapered, fenestrated, double stopcock	AR-3373-0302
Conical obturator w/ handle	AR-3375-0301
Sharp trocar w/ handle	AR-3375-0302
Conical obturator	AR-3375-0307
Switching stick, 3 mm × 38 cm	AR-3020-08

Inflow/Outflow Sheath for 3 mm Arthroscopes

I/O sheath set for Synergy 3 mm scope; part A & B	AR-3370-3003
Inner sheath for inflow/outflow Synergy 3 mm scope sheath	AR-3370-3003A
Outer sheath for inflow/outflow Synergy 3 mm scope sheath	AR-3370-3003B
Conical obturator	AR-3375-3005
Sharp trocar for inflow/outflow sheath	AR-3375-3006

2.7 mm Arthroscope

Arthroscope, 30°, 2.7 mm × 72 mm	AR-3350-2730
Arthroscope, 70°, 2.7 mm × 75 mm	AR-3350-2770
C-mount arthroscope, 30°, 2.7 mm × 72 mm	AR-3355-2730
C-mount arthroscope, 70°, 2.7 mm × 75 mm	AR-3355-2770

Sheaths and Obturators for 2.7 mm Arthroscopes

Sheath for 2.7 mm scopes, single stopcock	AR-3372-2701
Sheath for 2.7 mm scopes, double stopcock	AR-3372-2702
Conical obturator w/ handle	AR-3375-2701
Sharp trocar w/ handle	AR-3375-2702
Conical obturator	AR-3375-2707
Switching stick, 3.3 mm × 25 cm	AR-3020-09

2.4 mm Arthroscope

Arthroscope, 30°, 2.4 mm × 72 mm	AR-3350-2430
Arthroscope, 70°, 2.4 mm × 75 mm	AR-3350-2470
C-mount arthroscope, 30°, 2.4 mm × 72 mm	AR-3355-2430
C-mount arthroscope, 70°, 2.4 mm × 75 mm	AR-3355-2470

Sheaths and Obturators for 2.4 mm Arthroscopes

Sheath for 2.4 mm scopes, single stopcock	AR-3372-2401
Sheath for 2.4 mm scopes, double stopcock	AR-3372-2402
Conical obturator w/ handle	AR-3375-2401
Sharp trocar w/ handle	AR-3375-2402
Conical obturator	AR-3375-2407

1.9 mm Arthroscope


Arthroscope, 30°, 1.9 mm × 58 mm	AR-3350-1930
C-mount arthroscope, 30°, 1.9 mm × 58 mm	AR-3355-1930

Sheaths and Obturators for 1.9 mm Arthroscopes

Sheath for 1.9 mm scopes, single stopcock	AR-3372-1901
Sheath for 1.9 mm scopes, double stopcock	AR-3372-1902
Conical obturator w/ handle	AR-3375-1901
Sharp trocar w/ handle	AR-3375-1902
Conical obturator	AR-3375-1907

Centerline™ Arthroscope

Centerline carpal tunnel scope, 30°, 2.9 mm	AR-3350-2930
---	--------------



3.5 mm Hip Arthroscopes HD

Hip arthroscope, 30°, 3.5 mm × 202 mm	AR-3350-4030H
Hip arthroscope, 70°, 3.5 mm × 204 mm (a)	AR-3350-4070H
C-mount hip arthroscope, 30°, 3.5 mm × 202 mm	AR-3355-4030H
C-mount hip arthroscope, 70°, 3.5 mm × 204 mm (b)	AR-3355-4070H


3.5 mm Hip Sheaths

Sheath for 3.5 mm hip scopes, tapered, fenestrated, single stopcock	AR-3373-4001H
Sheath for 3.5 mm hip scopes, tapered, fenestrated, double stopcock	AR-3373-4002H
Conical obturator w/ handle	AR-3375-4001H
Sharp trocar w/ handle	AR-3375-4002H
Cannulated obturator w/ handle	AR-3375-4003H
Conical obturator	AR-3375-4007H
Switching stick, 3.5 mm × 38 cm	AR-3020-03
Switching stick, 3.5 mm × 46 cm	AR-3020-04
Switching stick, cannulated, 3.5 mm × 38 cm	AR-6520-16

Arthroscope Accessories

C-mount wrench	AR-3355-0001
----------------	--------------

Laparoscopes



3 mm Laparoscopes

Laparoscope, 0°, 3 mm × 207 mm	AR-3351-3000
Laparoscope, 30°, 3 mm × 207 mm	AR-3351-3030

5.5 mm Laparoscopes 4K

4K laparoscope, 0°, 5.5 mm × 300 mm (c)	AR-3351-5500
4K laparoscope, 30°, 5.5 mm × 300 mm	AR-3351-5530
4K laparoscope, 45°, 5.5 mm × 302 mm	AR-3351-5545
4K laparoscope, ext. length, 0°, 5.5 mm × 424 mm	AR-3351-5500L
4K laparoscope, ext. length, 30°, 5.5 mm × 424 mm	AR-3351-5530L
4K laparoscope, ext. length, 45°, 5.5 mm × 425 mm	AR-3351-5545L


5.5 mm High-Magnification Laparoscopes 4K

4K laparoscope, hi-mag, 0°, 5.5 mm × 300 mm	AR-3352-5500
4K laparoscope, hi-mag, 30°, 5.5 mm × 300 mm	AR-3352-5530
4K laparoscope, hi-mag, 45°, 5.5 mm × 302 mm	AR-3352-5545
4K laparoscope, hi-mag, NIR, 0°, 5.5 mm × 300 mm	AR-3352-5501
4K laparoscope, hi-mag, NIR, 30°, 5.5 mm × 300 mm	AR-3352-5531
4K laparoscope, hi-mag, NIR, 45°, 5.5 mm × 302 mm	AR-3352-5546
4K laparoscope, hi-mag, extended length, 0°, 5.5 mm × 424 mm	AR-3352-5500L
4K laparoscope, hi-mag, extended length, 30°, 5.5 mm × 424 mm	AR-3352-5530L
4K laparoscope, hi-mag, extended length, 45°, 5.5 mm × 425 mm	AR-3352-5545L


10 mm Laparoscopes 4K

4K laparoscope, 0°, 10 mm × 330 mm	AR-3351-1000
4K laparoscope, 30°, 10 mm × 330 mm (d)	AR-3351-1030
4K laparoscope, 45°, 10 mm × 333 mm	AR-3351-1045
4K laparoscope, 0°, 10 mm × 330 mm, NIR	AR-3351-1001
4K laparoscope, 30°, 10 mm × 330 mm, NIR	AR-3351-1031
4K laparoscope, 45°, 10 mm × 333 mm, NIR	AR-3351-1046
4K laparoscope, extended length, 0°, 10 mm × 455 mm	AR-3351-1000L
4K laparoscope, extended length, 30°, 10 mm × 455 mm	AR-3351-1030L
4K laparoscope, extended length, 45°, 10 mm × 459 mm	AR-3351-1045L

HD Endoscopes



HD endoscope, 0°, 4 mm × 174 mm (a)	AR-3000AW
HD endoscope, 30°, 4 mm × 176 mm	AR-3030AW
HD endoscope, 45°, 4 mm × 176 mm (b)	AR-3045AW
HD endoscope, 70°, 4 mm × 177 mm	AR-3070AW
HD endoscope, 30°, 4 mm × 176 mm, reverse light post	AR-3030AWR
HD endoscope, 45°, 4 mm × 176 mm, reverse light post	AR-3045AWR
HD endoscope, 70°, 4 mm × 177 mm, reverse light post (c)	AR-3070AWR



Open surgery Exoscope, 4K open scope 0° NIR, (a)	AR-3352-0050
Synergy Vision camera holder	AR-3210-0013


Scope and Sheath Accessories

Scope paste (rubbing compound/diamond paste)	8101001000
Stopcock assembly nut	8300018997
Stopcock handle (original design)	8300019001
Lubricant for stopcocks, Synergy	8300022404
Stopcock nut for handle	8300028607
Stopcock handle (new design)	8300028608
Replacement O-ring Kit for Synergy C-mount 5 per kit (black)	8300020859
Replacement O-ring Kit for Synergy stopcock 5 per kit (orange)	8300020860
Replacement O-ring Kit for bridge 5 per kit (clear)	8300026707
Black O-ring for metal cannula (5/pack)	8300031950
Replacement O-ring for outflow sheath (5/pack)	8300037687
Wolf light post adapter for scope	8300042617
Storz light post adapter for scope	89542801

Competitive Sheaths


Sheath, continuous flow, 2 stopcocks, rotatable, for Storz	AR-3300
Obturator, blunt, for AR-3300 sheath	AR-3301
Sheath, hi-flow, 2 stopcocks, rotatable, for Storz, tapered/fen	AR-3320-02
Sheath, hi-flow, 1 stopcocks, rotatable, for Storz, tapered/fen	AR-3320-01
Obturator, blunt, for AR-3320 series Storz sheath	AR-3321
Obturator, cannulated, for AR-3320 series Storz sheath	AR-3323
Trocar, for AR-3320 series Storz sheath	AR-3322
Obturator, conical, for AR-3320 series Storz sheath	AR-3324
Sheath, tap/fen, 2 stopcock, 4 mm Stryker	202-247-092
Trocar f/ 4 mm Stryker scope sheath	202-040-092
Obturator, conical f/ 4 mm Stryker scope sheath	202-042-092
Sheath, 2 stopcocks, 4 mm, 30 for Dyonics	203-247-106
Obturator, conical f/ 4 mm Dyonics scope sheath	203-042-106
Trocar f/ 4 mm Dyonics scope sheath	203-040-106

Instrument Cases



Video instrument tray	AR-3101C
Scope system case	AR-3102C
Tray for laparoscopes	AR-3103C
Tray for extended-length laparoscopes	AR-3104C
Midsize video instrument case (d)	AR-3105C
Tray for camera head	AR-3106C
Hip scope set instrument case	AR-3200C
Arthroscope case, for 2 scopes (e)	AR-3100

Bridges



- › Maintain portals while moving the arthroscope
- › Connect inflow or outflow tubing
- › J-lock–type connector
- › Single and double stopcock options

Bridge for Pano Scopes	AR-3380-4201 AR-3380-4202
Bridges for Sealed Cannulas for Pano	AR-3380-0421 AR-3380-0422

Bridges for 4 mm Arthroscopes

Bridge for 4 mm scope, single stopcock	AR-3380-4001
Bridge for 4 mm scope, double stopcock (b)	AR-3380-4002

Bridges for Sealed Cannulas

Sealed bridge for 4 mm scope, single stopcock (a)	AR-3380-0401
Sealed bridge for 4 mm scope, double stopcock	AR-3380-0402

Bridges for 3 mm Arthroscopes

Bridge for 3 mm scope, single stopcock	AR-3380-3001
Bridge for 3 mm scope, double stopcock	AR-3380-3002

Metal Cannulas and Sets




- › Reusable metal cannulas
- › Maintain portals while moving the arthroscope
- › Connect inflow or outflow tubing
- › Fenestrated and nonfenestrated options
- › J-lock–type connector
- › Threaded portion can lock into certain shaver blades
- › Sealed and unsealed versions

Cannulas and Obturators for Sealed Bridges

Cannula, metal, for bridge with dam, 4.5 mm	AR-3415-4.5
Cannula, metal, for bridge with dam, 5 mm	AR-3415-5.0
Cannula, metal, for bridge with dam (a) , 5.5 mm	AR-3415-5.5
Accessory, dam for bridge	AR-3491-05
Obturator, blunt, for bridge with dam, 4.5 mm	AR-3416-4.5
Obturator, blunt, for bridge with dam, 5 mm	AR-3416-5.0
Obturator, blunt, for bridge with dam, 5.5 mm	AR-3416-5.5

Cannulas and Obturators for Sealed Bridges for Pano™ Scope

Cannula, metal, 4.75 mm, for bridge w/ dam, Pano scope	AR-3515-4.75
Cannula, metal, 5.0 mm, for bridge w/ dam, Pano scope	AR-3515-5.0
Cannula, metal, 5.5 mm, for bridge w/ dam, Pano scope	AR-3515-5.5
Obturator, blunt for 4.75 mm metal cannulas w/ dam, Pano scope	AR-3516-4.75
Obturator, blunt for 5.0 mm metal cannulas w/ dam, Pano scope	AR-3416-5.0
Obturator, blunt for 5.5 mm metal cannulas w/ dam, Pano scope	AR-3416-5.5



Cannula Sets, Nonfenestrated (AR-3420-X.XS*) (A)

Product Description	Qty	Item Number
Metal cannula (a)	2	AR-3420-X.X
Blunt obturator (b)	1	AR-3421-X.X
Trocar (c)	1	AR-3422-X.X
Inflow cannula adapter (d)	1	AR-3035
Inflow cannula adapter, w/ stopcock (e)	1	AR-3035L

Cannula Sets, Fenestrated (AR-3425-X.XS) (B)

Product Description	Qty	Item Number
Metal cannula, fenestrated (f)	2	AR-3425-X.X
Blunt obturator (g)	1	AR-3421-X.X
Trocar (h)	1	AR-3422-X.X
Inflow cannula adapter (i)	1	AR-3035
Inflow cannula adapter, w/ stopcock (j)	1	AR-3035L

*X.X = 3 / 4.5 / 4.75 / 5 / 5.5 / 6.5

Cannulas and Obturators

Metal cannula, 3 mm	AR-3420-3.0
Metal cannula, 4.5 mm	AR-3420-4.5
Metal cannula, 4.75 mm	AR-3420-4.75
Metal cannula, 5 mm	AR-3420-5.0
Metal cannula, 5.5 mm	AR-3420-5.5
Metal cannula, 6.5 mm	AR-3420-6.5
Metal cannula, fenestrated, 3 mm	AR-3425-3.0
Metal cannula, fenestrated, 4.5 mm	AR-3425-4.5
Metal cannula, fenestrated, 4.75 mm	AR-3425-4.75
Metal cannula, fenestrated, 5 mm	AR-3425-5.0
Metal cannula, fenestrated, 5.5 mm	AR-3425-5.5
Metal cannula, fenestrated, 6.5 mm	AR-3425-6.5
Blunt obturator, 3 mm	AR-3421-3.0
Blunt obturator, 4.5 mm	AR-3421-4.5
Blunt obturator, 4.75 mm	AR-3421-4.75
Blunt obturator, 5 mm	AR-3421-5.0
Blunt obturator, 5.5 mm	AR-3421-5.5
Blunt obturator, 6.5 mm	AR-3421-6.5
Trocar, 3 mm	AR-3422-3.0
Trocar, 4.5 mm	AR-3422-4.5
Trocar, 4.75 mm	AR-3422-4.75
Trocar, 5 mm	AR-3422-5.0
Trocar, 5.5 mm	AR-3422-5.5
Trocar, 6.6 mm	AR-3422-6.5

Cannulas and Obturators for Pano Scope

Cannula, metal, 4.75 mm, Pano scope	AR-3420-4.75
Cannula, metal, 5 mm, Pano scope	AR-3420-5.0
Cannula, metal, 5.5 mm, Pano scope	AR-3420-5.5
Cannula, metal, 6.5 mm, Pano scope	AR-3420-6.5
Cannula, metal, 4.75 mm, fenestrated, Pano scope	AR-3425-4.75
Cannula, metal, 5 mm, fenestrated, Pano scope	AR-3425-5.0
Cannula, metal, 5.5 mm, fenestrated, Pano scope	AR-3425-5.5
Cannula, metal, 6.5 mm, fenestrated, Pano scope	AR-3425-6.5
Obturator, blunt for metal cannula, 4.75 mm, fenestrated, Pano scope	AR-3421-4.75
Obturator, blunt for metal cannula, 5 mm, fenestrated, Pano scope	AR-3421-5.0
Obturator, blunt for metal cannula, 5.5 mm, fenestrated, Pano scope	AR-3421-5.5
Obturator, blunt for metal cannula, 6.5 mm, fenestrated, Pano scope	AR-3421-6.5

Metal Outflow Cannulas and Accessories

Metal cannula, 3.2 mm w/ stopcock	SPK03299
Metal cannula, 3.2 mm w/o stopcock	SPK03296
Obturator f/ 3.2 mm metal cannula	SPK03297
Trocar f/ 3.2 mm metal cannula	SPK03298

Hip Bridges

Bridge for 3.5 mm hip arthroscope, single stopcock	AR-3380-4001H
Bridge for 3.5 mm hip arthroscope, double stopcock	AR-3380-4002H
Bridge for 3.5 mm hip arthroscope, quick connect, single stopcock	AR-3380-4003H
Bridge for 3.5 mm hip arthroscope, quick connect, double stopcock	AR-3380-4004H
FlushFit atraumatic cannula, quick-connect bridge, Pano™ scope	AR-3380-4056H

Hip Cannulas and Obturators for Quick-Connect Hip Bridges

Cannula, metal, hip, 4.5 mm	AR-3420H-4.5
Cannula, metal, hip, 5 mm	AR-3420H-5.0
Cannula, metal, hip, 5.5 mm	AR-3420H-5.5
Cannula, metal, hip, 6.5 mm	AR-3420H-6.5
Obturator, cannulated, hip, 4.5 mm	AR-3421H-4.5
Obturator, cannulated, hip, 5 mm	AR-3421H-5.0
Obturator, cannulated, hip, 5.5 mm	AR-3421H-5.5
Obturator, cannulated, hip, 6.5 mm	AR-3421H-6.5

Video Cart




To recapture valuable OR real estate, the Synergy video cart features a lighter and sleeker design with a smaller footprint than its predecessor. The cart is shipped with modular shelves preinstalled for recommended Arthrex component placement, so setup is quick and easy. The power cords are prewired and cable management can be accomplished on either side of the cart. Large wheels and handles on each side of the cart ensure ease of movement throughout the OR.

The boom arm on the main monitor allows 20 in of reach and 12 in of vertical adjustment. Arm tension is adjustable to hold 27 to 32 in monitors, and an easy-access handle and locking position facilitate cart movement. Mounted in the middle rear of the top shelf, the secondary arm is easily accessible by OR staff.

The base cart includes five shelves, a foot pedal holder, a secondary arm, and a main monitor boom arm.

- › Multiple configurations
- › Adjustable shelves
- › Rotatable touch panel arm

Video Cart and Monitor Stand



Multidirectional Boom Monitor Mount for 32 in and 17 in

IV Bag Hangers on Either Side (optional)

Fully Adjustable Shelves

Footswitch Holder

Smaller Footprint Than Previous Model

Touch Panel Arm (swings left or right)

Power Switch

Multisource Cable Connections in Rear

Recessed Top Shelf to Allow Loading Pump Tubing

CO₂ Bottle Holders (optional)

- > Multiple configurations
- > Adjustable shelves
- > 20 in extended boom arm
- > Rotatable touch panel arm



Accessories	
Video cart w/ boom, 120 V (a)	AR-3260-0001
Video cart w/ boom, 230 V	AR-3260-0002
Component shelf for video cart	AR-3260-1000
Drawer for video cart	AR-3260-1001
D & E size tank holder for video cart	AR-3260-1004
Footswitch bracket, AR-S8316C for video cart	AR-3260-1005
Footswitch bracket, AR-8310sp, and AR-8315w, for video cart	AR-3260-1006
Footswitch bracket, AR-8315c, for video cart	AR-3260-1007
Footswitch bracket, AR-6483 and AR-9800, for video cart	AR-3260-1008
Footswitch bracket, AR-8310, for video cart	AR-3260-1009
Dual-monitor assembly for video cart	AR-3260-1012
US bed rail for video cart	AR-3260-1014
Pump pullout shelf for video cart	AR-3260-1015
Cart handle	AR-3260-1017
IV bag hanger, 4 bags, for video cart	AR-3260-1023
Dual-yoke CO ₂ switching valve w/ hardware kit	AR-3260-1026
12G SDI, Isolated	NBB75DFIX
Connector, HDMI 2.0, F-F	7-004268-30
Splitter, C14 to dual C13	7-004269-30
Splitter, C14 to dual NEMA 5-15R	7-004270-30
Spare Parts	
Monitor handle	AR-3260-1018
Rear door, right, for video cart	AR-3260-1019
Rear door, left, for video cart	AR-3260-1020
Monitor boom arm for video cart	AR-3260-1021
Secondary arm for video cart	AR-3260-1022
Casters for video cart	AR-3260-1024
Casters Kit, 4, for video cart	AR-3260-1025
Patch panel	AR-3260-1029
3 pack paint, light gray, dark gray, black	AR-3260-1036
Hardware Kit for dual-yoke switching valve	AR-3260-1027
Housing, transformer	AR-3260-1037
Monitor dock	AR-3260-1038
Power cord, US	AR-3260-1040
Mushrooms (qty 2)	AR-3260-1055
Hinges (qty 2)	AR-3260-1057
Grommets, assorted	AR-3260-1056
Door magnet (qty 1)	AR-3260-1053
Channel nut screw (qty 6)	AR-3260-1054
Transformer, 120 V	AR-3260-1050

Secondary Monitor Stand

- > Small footprint
- > Adjustable monitor height
- > Large casters for ease of mobility
- > Locking casters
- > Dimensions 22 in (w) × 57.5 in (h) × 23 in (l)

Secondary monitor stand (a)	TP.5121.999
Monitor stand shelf w/ hardware	NT.4100.991
Monitor stand drawer w/ hardware	NT.4104.903
Monitor stand cable channel	NT.4115.806



Monitor Cart




Compact-Cart, 50U, 55" Monitor TP.5112.999

Printer and Recorder



Sony printer (b)	UP-DR80MD
Sony printer w/ wireless print system (printer and wireless combined)	UP-DR80MD/WU
Sony wireless print system accessory	UPA-WU10
Media pack, letter-size paper and ribbon	UPC-R81MD
Sony 4K video recorder	HVO-4000MT
Sony media tray, part 1 of 1	A1643186N
Sony media tray, part 2 of 2	A1643192H
Sony paper holder, left side (blue)	A5002529A
Sony paper holder, right side (pink)	A5002530A



Nano Arthroscopy

28	Nano Arthroscopy
29	NanoNeedle™ Scope Visualization System
33	NanoNeedle™ Scope Sheath and Access Kits
34	NanoScopic™ Release System
35	NanoNeedle™ Scope FiberTak® Anchor Integrated Sheath Guide
36	Nano Arthroscopy 2 mm Single-Use and Reusable Instrumentation
38	2.8 mm Nano Suction Punches
39	Nano Bone Prep Instrumentation
40	NanoResection™ Device Options
41	Fluid Management
42	NanoSuite
43	NanoSuite: Customizable Components

Nano Arthroscopy



The state-of-the-art NanoNeedle™ scope operative arthroscopy system uses high-quality, chip-on-tip image sensor technology to provide surgeons with a needle-sized, single-use camera system. The portability of the imaging system allows the surgeon to perform minimally invasive arthroscopy in the operating room using the NanoScope™ system as the main camera or as an adjunct to a traditional camera in a treatment room or in the physician's office.

A pioneer in Nano arthroscopy, Arthrex has designed the next generation of needle scope visualization systems. The NanoNeedle scope is a significant improvement in surgical imaging ergonomics and visualization quality.

The NanoNeedle scope provides an alternative to MRI imaging and second-look arthroscopy and offers precise image-guided visualization of injections.

The NanoNeedle kit and system seamlessly connects with the Synergy Vision™ or Synergy^D™ 4K system for full OR and EHR image management integration.

NanoNeedle™ Scope Visualization System



The NanoNeedle visualization system is a medical-grade, 3-in-1, chip-on-tip disposable camera system. With an intuitive tablet control unit, the system features the latest technologies in 1 mm imaging sensors, LED lighting, image management, and OR integration. A network-based system allows for bidirectional communication to your facility EHR, PACS, Synergy Surgeon™ app, and SurgeonVault® system.

NanoScope™ Console Specifications

- › Medical-grade camera control unit and camera card edge
- › 13 in, 3-in-1 camera control unit
- › Network capabilities to connect to facility's health record, PACS, and Synergy Surgeon™ app
- › Built-in microphone for video dictation without an external microphone
- › HDMI output to extend the video signal to in-room displays and integration systems
- › Printing capabilities
- › Adjustable brightness

NanoNeedle Scope 1.0

- › 400 × 400 resolution with 120° field of view
- › 5-50 mm depth of field
- › 125 mm, 180 mm, and 250 mm lengths
- › Compatible with NanoScope console and Synergy Vision™ console

NanoNeedle Scope 2.0

- › 720 × 720 resolution with 120° field of view
- › 5-50 mm depth of field
- › 125 mm, 180 mm, and 250 mm lengths
- › Compatible with Synergy Vision console


Economic Impact

Nano arthroscopy is a cost-effective alternative to standard arthroscopy. A single-use camera component eliminates procedure delays due to equipment cleaning, processing, and sterilization without costly maintenance, repairs, or upgrades related to traditional video stacks. The minimally invasive approach and unlimited access to joint spaces make the NanoScope system the instrument of choice for less invasive arthroscopic procedures.

Tablet Control Unit Components

- › 13 in HD monitor
- › Handpiece connector
- › Microphone
- › Ethernet, USB, and HDMI ports





Synergy Vision™ console	AR-3200-0026
Synergy Vision connect console	AR-3200-0027
NanoNeedle™ scope, 250 mm, single use	AR-3210-0045
NanoNeedle scope, 180 mm, single use	AR-3210-0044
NanoNeedle scope, 125 mm, single use	AR-3210-0043
NanoNeedle scope 2.0, 250 mm, single use (a)	AR-3210-0072
NanoNeedle scope 2.0, 180 mm, single use (b)	AR-3210-0071
NanoNeedle scope 2.0, 125 mm, single use (c)	AR-3210-0070
NanoScope™ tablet control unit	AR-3200-0030
Mobile cart (d)	AR-3502-CRT
NanoScope console battery supply replacement	150-0012-00-A
Monitor secondary stand	TP.5121.999
Monitor shelf w/ hardware	NT.4100.991
27 in LG display for secondary cart	27HQ710S
NanoScope probe	AR-10100N
Graduated black probe, reusable, 2.5 mm tip	AR-5007



NanoNeedle™ Scope Visualization System

The NanoNeedle visualization system is a medical-grade, 3-in-1, chip-on-tip disposable camera system. With an intuitive tablet control unit, the system features the latest technologies in 1 mm imaging sensors, LED lighting,

image management, and OR integration. A network-based system allows for bidirectional communication to your facility EHR, PACS, Synergy Surgeon™ app, and SurgeonVault® system.




NanoScope™ Console Specifications

- › Medical-grade camera control unit and camera card edge
- › 13 in, 3-in-1 camera control unit
- › Network capabilities to connect to facility's EHR, PACS, and Synergy Surgeon app
- › Built-in microphone for video dictation
- › HDMI output to extend the video signal to in-room displays and integration systems
- › Printing capabilities and brightness adjustment

Synergy Vision™ Console Specifications

- › Medical-grade camera control unit and camera card edge
- › 32 in, HDR surgical display
- › Network capabilities to connect to facility's EHR, PACS, and Synergy Surgeon app

NanoNeedle™ Scope Sheath and Access Kits



Introducing the NanoNeedle scope sheath and access kits, available in both low- and high-flow sheath options and a variety of lengths based on the area of use and applications. These comprehensive kits integrate seamlessly with FiberTak® anchors and can be ordered separately from the NanoNeedle scope.

NanoNeedle scope diagnostic sheath kit, 125 mm (bent and straight, 2.4 mm outer diameter)	AR-3210-0050
NanoNeedle scope high-flow operative sheath kit, 125 mm (bent and straight, 3.34 mm outer diameter) (a)	AR-3210-0051
NanoNeedle scope high-flow operative sheath kit (w/ crown tip and fenestration), 125 mm	AR-3210-0052
NanoNeedle scope sheath kit (high-flow access kit w/ spinal needle and guidewire), 125 mm	AR-3210-0053
NanoNeedle scope diagnostic sheath kit, 125 mm (bent and straight, 2.4 mm outer diameter)	AR-3210-0054
NanoNeedle scope high-flow operative sheath kit, 180 mm	AR-3210-0055
NanoNeedle scope sheath kit (diagnostic hip access kit w/ spinal needle and guidewire), 180 mm	AR-3210-0056
NanoNeedle scope sheath kit (high-flow hip access kit w/ spinal needle and guidewire), 180 mm	AR-3210-0057
NanoNeedle scope high-flow operative sheath kit (w/ crown tip and fenestration), 180 mm	AR-3210-0058
NanoNeedle scope working cannula, 1.5 cm	AR-3210-0059
NanoNeedle scope working cannula, 5 cm	AR-3210-0063
NanoNeedle diagnostic sheath kit, 250 mm, 2.4 mm outer diameter	AR-3210-0064
NanoNeedle high flow operative sheath kit, 250 mm, 3.34 mm outer diameter	AR-3210-0065
Portal skid	AR-4505
NanoGrip	AR-1090NH

NanoScopic™ Release System



A simplified, all-in-one system, the NanoScopic carpal tunnel release system streamlines ECTR procedures. Designed for precision and efficiency, this disposable system uses a retrograde cutting approach with a single-handed pull-blade technique to deliver a controlled carpal tunnel release, while its ergonomic design enhances comfort, feel, and overall instrument performance during the procedure.

From the operating room to the procedure room, experience the latest innovation in endoscopic carpal tunnel release with the NanoScopic carpal tunnel release system.

This revolutionary system promotes a faster procedural experience with ease, convenience, and mobility. Key highlights of the system include:

- > No trays or sterilization
- > No light guide
- > No coupler
- > No heat
- > No fogging


Increase your case volume with the mobility of the NanoScopic carpal tunnel release system.

NanoScopic™ carpal tunnel release kit (a)


AR-8850DS



NanoNeedle™ Scope FiberTak® Anchor Integrated Sheath Guide



DEX FiberTak handle, 125 mm, high-flow, crown-tip, straight w/ AR-8991 DX Knotless FiberTak anchor




Shoulder FiberTak handle, 125 mm, high-flow, crown-tip, straight w/ AR-3636 FiberTak anchor

Explore the first fully specialized Arthrex sheath to offer direct visualization of anchor insertion and placement for use in arthroscopic procedures. Anchor insertion can be achieved using the Nano-integrated sheath guides.

Arthrex offers 3 specialized guides—sized and calibrated to ensure an accurate depth stop for FiberTak anchor insertion—to accommodate distal extremity, shoulder, and hip procedures.


DX Knotless FiberTak anchor	AR-8991
Shoulder FiberTak anchor	AR-3636
NanoNeedle scope high-flow operative sheath kit, (w/ crown tip and fenestration), 125 mm	AR-3210-0052
NanoNeedle scope high-flow operative sheath kit, (w/ crown tip and fenestration), 180 mm	AR-3210-0058
Nano 1.8 FiberTak handle, shoulder	AR-3210-0060
Nano 1.8 FiberTak handle, DEX	AR-3210-0061
Bender	AR-4585

Nano Arthroscopy 2 mm Single-Use and Reusable Instrumentation



Harnessing 20 years of engineering excellence in designing arthroscopic hand instrumentation, Arthrex has produced the next generation in tissue resection and extraction instruments that are sharp and strong enough to resect and remove meniscal tissue.

Additionally, Arthrex engineered the smallest-diameter arthroscopic instrumentation for meniscal resection at just 2 mm in diameter. Nano instrumentation provides greater access for treating the most challenging pathology without sacrificing resection efficiency.



Knee, Hip, and Shoulder Nano Instrumentation, 130 mm Shaft Length, Single-Use Sterile Pack and Reusable Pack

NanoGrasper, straight, 130 mm, disposable (single pack)	AR-10913D-1
NanoScissor, straight, 130 mm, disposable (single pack)	AR-10915D-1
NanoBiter, straight, 130 mm, disposable (single pack)	AR-10911D-1
NanoBiter, 15° up, 130 mm, disposable (single pack)	AR-10922D-1
NanoBiter straight tip, 130 mm straight shaft, bendable, w/ FlushPort (a)	AR-10911F-1
NanoGrasper straight tip, 130 mm straight shaft, bendable, w/ FlushPort (b)	AR-10913F-1
NanoRetriever tip, 130 mm straight shaft, bendable, w/ FlushPort (c)	AR-10914F-1
NanoScissor tip, 130 mm straight shaft, bendable, w/ FlushPort (d)	AR-10915F-1
Nano BirdBeak® retriever tip, 130 mm straight shaft, bendable, w/ FlushPort (e)	AR-10916F-1
Nano BirdBeak grasper tip, 130 mm straight shaft, bendable, w/ FlushPort (f)	AR-10917F-1
Bending tool	AR-10900

Small Joint Nano Instrumentation, 70 mm Shaft Length, Single-Use Sterile Pack and Reusable Pack

NanoBiter, straight, 70 mm, disposable (single pack)	AR-10901D-1
NanoBiter, 15° up, 70 mm, disposable (single pack)	AR-10902D-1
NanoGrasper, straight, 70 mm, disposable (single pack)	AR-10903D-1
NanoScissor, straight, 70 mm, disposable (single pack)	AR-10905D-1
NanoBiter straight tip, 70 mm straight shaft, bendable, w/ FlushPort (g)	AR-10901F-1
NanoGrasper straight tip, 70 mm straight shaft, bendable, w/ FlushPort (h)	AR-10903F-1
NanoRetriever tip, 70 mm straight shaft, bendable, w/ FlushPort (i)	AR-10904F-1
NanoScissor tip, 70 mm straight shaft, bendable, w/ FlushPort (j)	AR-10905F-1
Nano BirdBeak retriever tip, 70 mm straight shaft, bendable, w/ FlushPort (k)	AR-10906F-1
Nano BirdBeak grasper tip, 70 mm straight shaft, bendable, w/ FlushPort (l)	AR-10907F-1
Bending tool (m)	AR-10900

2.8 mm Nano Suction Punches

Ideal for Nano arthroscopy, the Nano suction punch allows surgeons to achieve optimal resection and aspiration. Available in 2 sizes with 3 different cutting window configurations, this device can rotate 360° for up-, side-, or down-window resection.

Reusable and autoclavable, the Nano suction punch can accommodate all joints from shoulder to wrist applications. Moreover, the Nano suction punch integrates with the GraftNet™ device to collect aspirated tissue, making it ideal for tissue biopsies, AutoCart™ procedures, and autograft biologic augmentation.

Nano suction punch, 70 mm flat tip (a)	AR-10920F-1
Nano suction punch, 70 mm scoop tip (b)	AR-10921F-1
Nano suction punch, 70 mm bullet tip (c)	AR-10922F-1
Nano suction punch, 130 mm flat tip (d)	AR-10930F-1
Nano suction punch, 130 mm scoop tip (e)	AR-10931F-1
Nano suction punch, 130 mm bullet tip (f)	AR-10932F-1



Flat Tip **(a, d)**



Scoop Tip **(b, e)**



Bullet Tip **(c, f)**

Nano Bone Prep Instrumentation



With the bone prep instrumentation line, surgeons have the instrumentation necessary to complete efficient and effective microfracture and small bony work. Nano arthroscopy, a line that will continue to expand, can help surgeons grow their minimally invasive practices.

Key Features and Benefits

- > **Bendability:** Contour with the Nano reusable bending tool
- > **Sharpness:** Single-use instruments will always be sharp and precise
- > **Size:** 2 lengths (7 cm and 10 cm) ergonomically designed for expanded applications


Technical Pearls

- > **Compatibility:** Low-profile tools pair with AutoCart™ surgical technique for chondrocyte transplantation
- > **Efficiency:** Continue with the Nano Sabre shaver and GraftNet™ device and combine with a BioCartilage® allograft to complete the AutoCart procedure

30° chondral pick, 7 cm (a)	AR-1091CP-70-1
30° chondral pick, 10 cm (b)	AR-1091CP-100-1
30° chondral pick, 7 cm (5-pack)	AR-1091CP-70
30° chondral pick, 10 cm (5-pack)	AR-1091CP-100
Cup curette, 7 cm (c)	AR-1091CC-70-1
Cup curette, 10 cm (d)	AR-1091CC-100-1
Cup curette, 7 cm (5-pack)	AR-1091CC-70
Cup curette, 10 cm (5-pack)	AR-1091CC-100
Ring curette, 7 cm (e)	AR-1091RC-70-1
Ring curette, 10 cm (f)	AR-1091RC-100-1
Elevator, 7 cm (g)	AR-1091E-70-1
Elevator, 10 cm (h)	AR-1091E-100-1



NanoResection™ Device Options



NanoResection small-hub shaver handpiece, blades, and burrs allow for efficient minimally invasive tissue resection and debridement. Low-profile tips facilitate safe introduction into most tight joint spaces without a limb holder.

The NanoResection small-hub shaver handpiece is a lightweight, high-speed, and high-torque offering within the Synergy^{Resection}™ system. When combined with small-diameter attachments optimized for Nano arthroscopy, this shaver system offers the most versatile resection devices available.

The Apollo^{RF} SJ50 probe is the first optimized aspirating probe for Nano arthroscopy. Featuring 360° edge control, an anatomically designed 50° curve, and a low default ablation setting, the probe ensures precise ablation and coagulation while minimizing effect on adjacent tissue.


Synergy ^{Resection} console (a)	AR-8305
NanoResection small-hub shaver handpiece (b)	AR-8330SJ
Nano Sabre shaver (c)	AR-9280NSR
Nano oval burr (d)	AR-9280NOT
Low-profile footswitch	AR-8310
Cordless gas pedal footswitch (e)	AR-8315W
Corded gas pedal footswitch	AR-8315C
Synergy ^{RF} console (f)	AR-9800
Synergy ^{RF} foot pedal (g)	AR-9800-F
Apollo ^{RF} SJ50 probe (h)	AR-9845

Fluid Management



Arthrex offers a comprehensive portfolio of powered resection and fluid management products that provides options for all procedure types and settings. Whether it's controlling distention during complex arthroscopic procedures or providing simple gravity flow during office procedures, Arthrex is dedicated to offering unique solutions designed for safe and reliable joint distention, tissue resection, and small-profile visualization with Nano arthroscopy in the OR or treatment room.

DualWave™ arthroscopy pump (a)	AR-6480
Main pump tubing	AR-6410
DualWave outflow tubing (b)	AR-6430
DualWave inflow and outflow tubing	AR-6436
Waste collection bag	AR-6431
ReDeuce™ pump tubing	AR-6411
ReDeuce patient tubing	AR-6421
DualWave outflow tube set w/ ReDeuce tubing system	AR-6435
Extension tubing	AR-6220
Continuous Wave™ 4 arthroscopy pump	AR-6485
Continuous Wave 4 inflow tubing	AR-6413
Gravity tubing, 2 spike	AR-6412
Gravity tubing, 4 spike	AR-6414
DualWave remote (d)	AR-6495
DualWave foot pedal	AR-6483
DualWave pump cart (c)	AR-6481
NanoScope™ console mounting arm	AR-5010-1500
Monitor stand	002309-00
Outflow power strip	002313-00
Tablet arm for Nano console mount	304855-00
Adjustable shelf for stand	304813-01
Nano Arthroscopy Percutaneous Insertion Kit (includes calibrated spinal needle-stylet and 18 ga guidewire for direct cannula insertion)	AR-1090PK-1
NanoNeedle™ scope working cannula, 1.5 cm	AR-3210-0059
NanoNeedle scope working cannula, 5 cm	AR-3210-0063





Elevate Your Patients' Arthroscopy Experience

The NanoSuite is a completely customizable procedure room that creates a true same-day surgery experience. Surgeons can offer patients a minimally invasive option by using local anesthesia, facilitating an intimate dialogue that creates trust and leads to better outcomes.^{1,2} Facilities can leverage this technology to drive new business and help gain and retain patients.

With the list of equipment on the following page, your facility can offer the next revolution in orthopedics to diagnose and treat your patients less invasively. Facilities that have incorporated or adopted Nano arthroscopy technology experience a direct impact on

patient satisfaction and reduced capital outlay needed to perform arthroscopic techniques. Patients have reported higher satisfaction scores,² returned to work and play sooner,¹ and overwhelmingly referred their friends and family for a Nano procedure.²

By offering an all-disposable method, facilities may experience lower operating costs, save on sterilization, and notice lower anesthesia costs. All things combined, the NanoSuite is the complete solution for facilities to differentiate themselves from the competition, while offering the cutting edge of medical innovation.

References

1. Colasanti CA, Mercer NP, Garcia JV, Kerkhoffs GMMJ, Kennedy JG. In-office needle arthroscopy for the treatment of anterior ankle impingement yields high patient satisfaction with high rates of return to work and sport. *Arthroscopy*. 2021;S0749-8063(21)00848-3. doi:10.1016/j.arthro.2021.09.016
2. Murawski CD, Kennedy JG. Anteromedial impingement in the ankle joint: outcomes following arthroscopy. *Am J Sports Med*. 2010;38(10):2017-2024. doi:10.1177/0363546510369335

NanoSuite

NanoSuite: Customizable Components




Surgeons can tailor their procedure rooms to meet their unique needs. They may select from a variety of components that offer the latest technology while promoting a positive patient experience:

- › Consoles
- › Chairs and Accessories
- › Nano Camera Handpieces
- › Carts and Accessories
- › Nano Arthroscopy Imaging System
- › Nano Sheath Kit
- › 2.7 mm Nano Cannula and Insertion Kits
- › 2 mm Nano Arthroscopy Single-Use Instrumentation
- › Arthroscopy Powered Resection Product and Accessories
- › Fluid Management Products and Accessories
- › Cell Therapy Products

Synergy Power™ System

46	Synergy Power™ System Overview
47	Synergy Power™ Dual Trigger Rotary Drill
48	Synergy Power™ Sagittal Saw
49	Drilling Attachments
50	Burring Attachments
50	Reaming Attachments
50	Sawing Attachments
51	Saw Blades
54	Sterilization Containers
55	Battery Packs and Chargers

Synergy Power™ System Overview



The Synergy Power system is a powerful, versatile, and reliable battery-powered system designed for a multitude of orthopedic applications, including sports, arthroplasty, trauma, and distal extremities procedures. With two handpieces—the dual trigger rotary drill and dedicated sagittal saw—the Synergy Power system provides surgeons and facilities with precision, power, and flexibility. The dual trigger rotary drill can drill, ream, saw, and burr while delivering the speed needed for heavy drilling procedures and the necessary for heavy-duty reaming applications. Designed for efficiency, the

dedicated sagittal saw offers a robust and precise cutting experience with a proprietary open-hub design for easy visualization and confirmation of blade placement.

The Synergy Power system offers lithium-ion battery packs that are sterilizable. Designed for procedure customization, the sterilizable battery packs come in a small size for light-duty procedures and a large size for heavy-duty procedures that require maximum battery run time. Engineered to have a small footprint, a simple user interface, four charging bays, and quick charging times.

Synergy Power™ Dual Trigger Rotary Drill



Features and Benefits

- › Versatile performance for sports, arthroplasty, trauma, and distal extremity procedures
- › Modular handpiece to drill, ream, saw, and burr with power and precision
- › Ergonomic, lightweight aluminum body
- › Dual triggers for forward, reverse, and oscillation modes
- › IPX-rated 6, 8, and 9 for the demanding sterile processing environment
- › Innovative twist-collet mechanism to easily connect and remove attachments

Dual trigger rotary drill

AR-01-01-001

Synergy Power™ Sagittal Saw



























Features and Benefits



- > Powerful cutting performance across a range of surgical applications
- > Ergonomic, lightweight aluminum body
- > Proprietary open-hub design for clear visibility to the blade connection point
- > Indexing head with 8 cutting positions
- > IPX-rated 6, 8, and 9 for the demanding sterile processing environment

Sagittal saw	AR-01-01-002
--------------	--------------



Drilling Attachments

Drill Attachment	End Effect	Product Name	Description	Product Number
		Small AO drill	One rounded and one flat side	AR-01-02-004
		Hudson drill	Flowerpot-shaped locking mechanism	AR-01-02-006
		Hudson modified trinkle drill	Flowerpot-shaped locking mechanism; shaft with 3 flat surfaces	AR-01-02-005
		Trinkle drill	Round shank with ball indent	AR-01-02-007
		J/T-latch	Accepts T-latch and J-latch	AR-01-02-026
		1/4" drill chuck	Varying geometries between 0.5-7.4 mm	AR-01-02-018
		1/4" drill chuck, hybrid	Varying geometries between 0.3-7.4 mm	AR-01-02-020
		1/8" drill chuck, hybrid	Varying geometries between 0.3-4.0 mm	AR-01-02-021
		5/32" drill chuck	Varying geometries between 0.5-4.0 mm	AR-01-02-017
		Dual trigger wire collet, 0.6-1.8 mm	Accepts wire between 0.6-1.8 mm	AR-01-02-001
		Dual trigger pin collet, 2.0-3.2 mm	Accepts pin between 2.0-3.2 mm	AR-01-02-002
		Dual trigger pin collet, 3.0-4.2 mm	Accepts pin between 3.0-4.2 mm	AR-01-02-003

Burring Attachments


Bur Attachment	End Effect	Product Name	Description	Product Number
		High-speed bur	Accepts T-latch, J-latch, and round end effects with 2.35 mm diameter	AR-01-02-009

Reaming Attachments

Ream Attachment	End Effect	Product Name	Description	Product Number
		Hudson modified trinkle reamer	Flowerpot-shaped locking mechanism; shaft with 3 flat surfaces	AR-01-02-011
		Trinkle reamer	Round shank with ball indent	AR-01-02-013
		Large AO reamer	Round shank with horizontal groove	AR-01-02-014
		Hudson reamer	Flowerpot-shaped locking mechanism	AR-01-02-012
		DHS/DCS reamer	Shaft with 3 flat surfaces and vertical groove	AR-01-02-016
		1/4" reamer chuck	Varying geometries between 0.5-7.4 mm	AR-01-02-019

Sawing Attachments

Only Arthrex blades are compatible with the Synergy Power™ system's sawing attachments.

Saw Attachment	End Effect	Product Name	Description	Product Number
		Micro sagittal saw	Accepts a variety of micro sagittal saw blades	AR-01-02-024
		Reciprocating saw	Accepts a variety of reciprocating saw blades	AR-01-02-025

Saw Blades















Large Sagittal Blades with Dedicated Sagittal Saw Handpiece

Blade	Description	Length × Width × Cut Thickness/ Material Thickness	Part Number
	Large bone sagittal blade	90 mm × 12.5 mm × 0.89 mm	AR1-125-90-089-SAG
	Large bone sagittal blade	90 mm × 12.5 mm × 1.19 mm	AR1-125-90-119-SAG
	Large bone sagittal blade	90 mm × 12.5 mm × 1.27 mm	AR1-125-90-127-SAG
	Large bone sagittal blade	90 mm × 12.5 mm × 1.37 mm	AR1-125-90-137-SAG
	Large bone sagittal blade	90 mm × 19 mm × 0.99 mm	AR1-19-90-099-SAG
	Large bone sagittal blade	90 mm × 19 mm × 1.19 mm	AR1-19-90-119-SAG
	Large bone sagittal blade	90 mm × 19 mm × 1.27 mm	AR1-19-90-127-SAG
	Large bone sagittal blade	90 mm × 19 mm × 1.37 mm	AR1-19-90-137-SAG
	Large bone sagittal blade	90 mm × 21 mm × 1.19 mm	AR1-21-90-119-SAG
	Large bone sagittal blade	90 mm × 21 mm × 1.27 mm	AR1-21-90-127-SAG
	Large bone sagittal blade	90 mm × 21 mm × 1.37 mm	AR1-21-90-137-SAG
	Large bone sagittal blade	90 mm × 25 mm × 0.89 mm	AR1-25-90-089-SAG
	Large bone sagittal blade	90 mm × 25 mm × 0.99 mm	AR1-25-90-099-SAG
	Large bone sagittal blade	90 mm × 25 mm × 1.19 mm	AR1-25-90-119-SAG
	Large bone sagittal blade	90 mm × 25 mm × 1.27 mm	AR1-25-90-127-SAG
	Large bone sagittal blade	90 mm × 25 mm × 1.37 mm	AR1-25-90-137-SAG

■ DePuy Synthes Attune ■ Stryker Triathlon ■ Zimmer Persona

Blade recommendations for specific procedures

Large Sagittal Blades with Dedicated Sagittal Saw Handpiece

Blade	Description	Length × Width × Cut Thickness/ Material Thickness	Part Number
	Large bone sagittal blade	90 mm × 25 mm × 1.47 mm	AR1-25-90-147-SAG
	Large bone sagittal blade	100 mm × 19 mm × 1.19 mm	AR1-19-100-119-SAG
	Large bone sagittal blade	100 mm × 19 mm × 1.27 mm	AR1-19-100-127-SAG
	Large bone sagittal blade	100 mm × 19 mm × 1.37 mm	AR1-19-100-137-SAG
	Large bone sagittal blade	100 mm × 25 mm × 1.27 mm	AR1-25-100-127-SAG
	Large bone sagittal blade	75 mm × 25 mm × 1.27 mm	AR1-25-75-127-SAG
	Large bone sagittal blade	70 mm × 12.5 mm × 0.99 mm/0.64 mm	AR1-125-70-099-SAG
	Large bone sagittal blade	70 mm × 18.5 mm × 0.91 mm/0.79 mm	AR1-185-70-091-SAG
	Large bone sagittal blade	70 mm × 25 mm × 0.94 mm/0.89 mm	AR1-25-70-094-SAG
	Large bone sagittal blade	70 mm × 25 mm × 1.47 mm/1.27 mm	AR1-25-70-147-SAG
	Large bone sagittal blade	64 mm × 35 mm × 0.89 mm/0.80 mm	AR1-35-64-089-SAG
	Large bone sagittal blade	60 mm × 34.5 mm × 0.79 mm	AR1-345-60-079-SAG
	Large bone sagittal blade	60 mm × 35 mm × 1.35 mm/1.19 mm	AR1-35-60-135-SAG
	Large bone sagittal blade	60 mm × 48.5 mm × 0.99 mm/0.64 mm	AR1-485-60-099-SAG
	Large bone sagittal blade	55 mm × 29 mm × 1.02 mm/0.64 mm	AR1-29-55-102-SAG
	Large bone sagittal blade	30 mm × 9 mm × 0.79 mm/0.64 mm	AR1-9-30-079-SAG

■ DePuy Synthes Attune ■ Stryker Triathlon ■ Zimmer Persona

Blade recommendations for specific procedures

Reciprocating Blades with Reciprocating Saw Attachment



Blade	Description	Length × Width × Cut Thickness/ Material Thickness	Part Number
	Reciprocating saw blade (sternum blade)	32.8 mm × 9.5 mm × 1.1 mm	AR-600-025S
	Reciprocating saw blade	77.5 mm × 11.2 mm × 1.17 mm/0.76 mm	AR1-112-775-117RCP
	Reciprocating saw blade	77.5 mm × 11 mm × 1.60 mm/1.24 mm	AR1-11-775-160RCP
	Reciprocating saw blade	70.8 mm × 12.5 mm × 1.12 mm/1.0 mm	AR1-125-708-112RCP
	Reciprocating saw blade	70 mm × 12.5 mm × 0.94 mm/0.64 mm	AR1-125-70-094RCP
	Reciprocating saw blade	60 mm × 6.3 mm × 1.04 mm/0.64 mm	AR1-63-60-104RCP

Micro Sagittal Blades with Micro Sagittal Saw Attachment



Blade	Description	Length × Width × Thickness	Part Number
	Micro sagittal saw blade	40 mm × 14 mm × 0.6 mm	AR-300-040S
	Micro sagittal saw blade	40 mm × 9.5 mm × 0.6 mm	AR-300-041S
	Micro sagittal saw blade EZ-90	25 mm × 9.4 mm × 0.7 mm	AR-300-042S
	Micro sagittal saw blade	25 mm × 9.4 mm × 0.6 mm	AR-300-043S
	Micro sagittal saw blade	25 mm × 5.5 mm × 0.6 mm	AR-300-044S
	Micro sagittal saw blade	16 mm × 5.5 mm × 0.6 mm	AR-300-045S
	Micro sagittal saw blade	18.5 mm × 10 mm × 0.55 mm	AR-300-407S
	Micro sagittal saw blade with stop	17 mm × 9.4 mm × 0.55 mm	AR-300-406S
	Micro sagittal saw blade, angled	19 mm × 10 mm × 0.60 mm	AR-300-450S

Sterilization Containers




The Synergy Power™ sterilization containers are engineered for durability and ease of use, providing reliable protection and readiness for all your procedural needs.

Specifically designed to enhance workflow efficiency in operating rooms, these trays offer grab-and-go versatility, ensuring they are ready to support a wide array of procedures immediately upon removal from storage.

The Synergy Power system offers 2 sterilization tray configurations, each designed to support the full range of Synergy Power handpieces and attachments. Whether your team prefers a streamlined setup or expanded instrument access, both trays deliver uncompromising compatibility and efficiency.

2-Handpiece Tray


Designed for efficiency, optimized for streamlined setups



2-Handpiece container bottom	AR-01-05-004
2-Handpiece insert	AR-01-05-005
2-Handpiece container lid	AR-01-05-006

3-Handpiece Tray

Expanded instrument capacity for maximum surgical flexibility



3-Handpiece container bottom	AR-01-05-007
3-Handpiece insert	AR-01-05-008
3-Handpiece container lid	AR-01-05-009

Battery Packs and Chargers



The Synergy Power™ system is powered by 13.2 V lithium-ion battery packs, enabling the handpieces to deliver optimal power for peak performance. The availability of both small and large battery packs gives facilities the flexibility to choose the best option for each procedure. Designed as an intuitive, space-saving solution, Synergy Power lithium-ion battery chargers boast smart charging capabilities with simple iconography for charge status at a glance.

Lithium-ion battery packs (a)

- › 13.2 V lithium-ion batteries deliver optimal power for peak performance
- › Sterilizable batteries offered in 2 sizes provide flexibility for varying caseloads
- › Sterilizable and aseptic* options allow facilities to choose their ideal battery solution

Lithium-ion battery chargers (b)

- › Smart chargers designed to provide a fast and efficient charge while taking up minimal space
- › Intuitive design for clear confirmation of battery charge level
- › Charge 4 batteries at a time

Lithium-ion battery, small ST	AR-01-03-BP
Lithium-ion battery, large ST	AR-01-03-LB
Battery charger for lithium-ion batteries	AR-01-03-001

*Product not yet available for sale

Synergy^{Resection™} System

58	Synergy ^{Resection™} Console
59	Handpieces
60	Footswitches
61	Large-Hub Options
63	Retractable Hood Burrs
63	PowerPick™ and PoweRasp™ Attachments
64	NanoResection™ Shaver Handpiece Options (7 cm)
64	NanoResection™ Instruments (11 cm length)
65	ShaverDrill™ Instruments

Synergy^{Resection}™ Console



This versatile Synergy^{Resection} system was designed for multiple orthopedic procedures and operates like 2 consoles in 1 chassis. It can simultaneously operate 2 accessory handpieces completely independent of each other. Operators can also use 1 or 2 foot pedals to control each of the accessory handpieces. Shaver handpieces offer 3 oscillating modes—standard, efficient, and aggressive—depending on the application or surgeon’s preference.

Integration of the Synergy^{Resection} console with a Synergy video management system enables the heads-up display feature and allows for easy visibility of the shaver system settings on the monitor.

When the Synergy^{Resection} console is integrated with the DualWave™ or Continuous Wave™ 4 arthroscopy pump, the pump detects when the shaver handpiece is active and automatically compensates for the loss of distention by increasing the pressure of the active shaver.

- › Multifunctional shaver system
- › Operate 2 handpieces simultaneously
- › Automatic detection of handpiece
- › Control via 1 or 2 foot pedals
- › 3 different oscillating modes: standard, efficient, and aggressive
- › Large and small hub: OSC 3000 rpm max
- › Large hub: FWD/REV 8000 rpm max
- › Small hub: FWD/REV 6200 rpm max

Synergy^{Resection} console

AR-8305

Handpieces



- Hand Control**
- > Multifunction buttons
 - > Change speed
 - > Change direction
 - > Turn on/off

(a)



- Hand Control (backhand)**
- > Multifunction buttons
 - > Change speed
 - > Change direction
 - > Turn on/off
 - > Suction opposite from button controls

(b)



- Foot Control**
- > Lightweight
 - > Control via footswitch

(c)




- NanoResection Small-Hub Shaver Handpiece**
- > Maximum control and dexterity
 - > Compact and lightweight
 - > Control via footswitch
 - > Used in Nano applications
 - > Accepts NanoResection and small-hub shaver blades and burrs

(d)



(e)

(f)



(g)



(h)

Arthrex offers unique double-sealed shaver handpieces that are robust and reliable with a light ergonomic feel, with or without hand-control options. These high-speed, high-torque handpieces feature a snap-lock collet that accepts standard, hip-length, and small joint blades and burrs, making it one of the most versatile handpieces available. With high-performance motors and constructed of cutting-edge materials, Synergy^{Resection}™ shaver handpieces are a great addition to the OR.

Hand-Control Shaver Handpieces

Featuring ergonomic multifunction buttons that place control of speed and direction at the surgeon's fingertips, the hand-control shaver handpieces are available in standard and backhand models.

Foot-Control Shaver Handpiece

With the same lightweight, high-performance features as our hand-operated handpieces, the footswitch shaver handpiece provides another option for control.

NanoResection™ Small-Hub Shaver Handpiece

This lightweight and ergonomic handpiece accepts small-hub attachments and provides highly precise, controlled resection.

Shaver handpiece, hand control (a)	AR-8332H
Shaver handpiece, hand control, backhand (b)	AR-8332RH
Shaver handpiece, footswitch control (c)	AR-8330F
NanoResection small-hub shaver handpiece (d)	AR-8330SJ
Replacement valve assembly for AR-8330F, AR-8330H, AR-8330RH, AR-8332H, and AR-8332RH handpieces (e)	AR-8330V
Replacement O-rings for AR-8330V valve assembly (pack of 10) (f)	AR-8330VO
Replacement O-rings for F-style soaker cap (pack of 5) (g)	AR-8330FO
Replacement O-rings for T-style soaker cap (pack of 5) (h)	AR-8330TO
Replacement valve assembly for AR-8330SJ	AR-8330JV
Replacement O-ring for AR-8330JV valve assembly (pack of 10)	AR-8330SJO

Footswitches



Low-Profile Footswitch

The low-profile footswitch minimizes user fatigue while maintaining control of all functions. Control direction and speed by pressing the appropriate pad.

Press the toggle dome switch to change control between the console's dual channels if 2 handpieces are plugged in. If 1 handpiece is plugged in, the toggle dome switch will toggle between the 3 oscillation modes.

Wireless and Corded Multifunction Footswitches

The wireless and corded multifunction footswitches feature user-selectable functionality with either standard footswitch or gas-pedal modes. All speed, direction, cruise, and toggle controls are accessible on either of these user-friendly accessories, and the active mode is displayed on the console.

Low-profile footswitch (e)	AR-8310
Synergy ^{Resection} ™ wireless footswitch (f)	AR-8315W
Corded gas-pedal footswitch (g)	AR-8315C












Large-Hub Options

Fits Large Shaver Handpieces

	2 (Ø mm)	3 (Ø mm)	3.5 (Ø mm)	3.8 (Ø mm)	4 (Ø mm)	4.2 (Ø mm)	5 (Ø mm)	5.5 (Ø mm)
 <p>Torpedo™ Blade Tapered tip and scissor-like cutting action make for easy access into joint space and rapid soft-tissue resection; outer oval window to limit tissue size and help reduce clogging</p>			AR-7350TD ¹		AR-7400TD ¹ AR-8400TD AR-8400CTD (curved)	AR-6420TD ² AR-6420CTD ³	AR-8500TD	
 <p>Excalibur Sharp-toothed blade makes this the most aggressive shaver for extensive soft-tissue debridement</p>				AR-8380EX	AR-8400EX AR-8400CEX (curved)	AR-6420EX ² AR-6420CEX ³	AR-8500EX	AR-8550EX
 <p>Bone Cutter Dual-purpose blade designed for fast and aggressive tissue debridement when in oscillation mode and for bone removal when run in the forward or reverse direction</p>				AR-8380BC	AR-8400BC AR-8400CBC (curved)	AR-6420BC ² AR-6420CBC ³	AR-8500BC	AR-8550BC
 <p>Double Cut Dual inner-cutting windows provide efficient soft-tissue resection; elongated outer window minimizes tissue size entering the blade, thus reducing clogging</p>				AR-8380DC	AR-8400DC AR-8400CDC (curved)			AR-8550DC
 <p>SabreTooth Wide-toothed blade for bulk soft-tissue debridement</p>			AR-8350ST	AR-8380ST	AR-8400ST AR-8400CST (curved)	AR-6420CST ³	AR-8500ST	AR-8550ST
 <p>Sabre Smooth teardrop-shaped outer window for general soft-tissue resection that leaves smooth anatomic edges</p>	AR-7200SR ¹	AR-7300SR ¹		AR-8380SR	AR-8400SR AR-8400CSR (curved)		AR-8500SR	AR-8550SR
 <p>Dissector Sharp inner-toothed blade with a teardrop-shaped outer window for less iatrogenic chondroplasties and soft-tissue removal</p>		AR-7300DS ¹	AR-8350DS AR-8350CDS (curved)	AR-8380DS	AR-8400DS AR-8400CDS (curved)	AR-6420DS ² AR-6420CDS ³	AR-8500DS	AR-8550DS
 <p>Slotted Whisker Double-slotted outer window enables precise tissue resection</p>				AR-8380SW				

Color bars reflect family color shown on the adapter



1. Small-joint length, large hub
2. Hip length, straight
3. Hip length, curved

		3 (Ø mm)	4 (Ø mm)	5 (Ø mm)	5.5 (Ø mm)	
Color bars reflect family color shown on the adapter	Diamond		Hip DiamondCut™ Burrs Designed for precise control during resection of femoroacetabular impingement (FAI) associated with pincer and cam deformities, the diamond-coated burr head enables smooth bone resection that helps control and polish		AR-6400RBD	AR-6500RBD
	2-flute		BLURR® Attachment A hybrid burr and blade for aggressive resection of bone and soft tissue			AR-8500BL
	6-flute 8-flute 12-flute		ClearCut™ Oval Burrs Oval burr with clear hood to aid in providing maximum visibility and protecting surrounding anatomy		AR-8500COE AR-8500COT	AR-8550COS AR-8550COE AR-8550COT
	8-flute		ClearCut Round Burrs Round burr with clear hood to aid in providing maximum visibility and protecting surrounding anatomy		AR-8400CRE	AR-8500CRE AR-8550CRE
	6-flute		ClearCut TR Burrs Angled, clear hood to retract and protect sensitive tissue while burring		AR-8400CTS	
	6-flute		ClearCut SLAP Burrs Clear hood fully encases the burr head for increased protection of surrounding tissue in SLAP lesion repairs of the labrum		AR-8400CSS	
	6-flute 8-flute 12-flute		FlushCut™ Oval Burrs Burr head is flush with the outer tube to create an even cutting surface, which helps to prevent troughs and scalloping in bone		AR-8500FOS AR-8500FOE AR-8500FOT	AR-8550FOS AR-8550FOE AR-8550FOT
	8-flute 10-flute 12-flute		Oval Burrs Designed for rapid, aggressive bone resection in notchplasty, distal clavicle, and subacromial decompression procedures		AR-7300OBT ¹	AR-8400OBE AR-8500OBE AR-8550OBT
	8-flute		Round Burrs A precise and aggressive burr designed for notchplasties, glenoid and acetabulum preparation, and osteophyte removal		AR-7300RBE ¹	AR-8400RBE AR-6400RBE ² AR-8500RBE AR-6500RBE ² AR-6500VBE ³ AR-8550RBE
	12-flute		Tapered Burrs Excellent burr for distal clavicle resection and for shaping bone and cartilage during notchplasty and acromioplasty procedures			AR-8500TBT
		Hoodless Sheath, O/T Burrs Accessory outer tube to fully expose the burr head for access to tight spaces and where maximum visibility is desired		AR-8400HSO oval/burr	AR-8500HSO oval/ tapered burr AR-8550HSO oval/burr	

1. Small-joint length, large hub
2. Hip length, straight
3. Hip length, curved



Retractable Hood Burrs


Fits Large-Hub Shaver Handpieces

		4 (Ø mm)	5 (Ø mm)	5.5 (Ø mm)
Color bars reflect family color shown on the adapter	8-flute	 <p>Retractable FlushCut™ Oval Burrs Designed to fully expose the burr head intraoperatively with a single-handed flick of a lever; ideal for distal clavicle resection where maximal burr head exposure is needed and for acromioplasties, where the hood can guard sensitive tissue from the burr head</p>	AR-8500RFOE	AR-8550RFOE
	8-flute	 <p>Retractable Hip-Length Round Burrs Uniquely designed burr with the ability to fully expose the burr head intraoperatively with a single-handed flick of a lever; ideal for cam resections where maximal burr-head exposure is needed, and for acetabuloplasties, where the hood can guard sensitive tissues</p>	AR-6400RRBE	AR-6500RRBE AR-6550RRBE








PowerPick™ and PowerRasp™ Attachments

Fits Large Shaver Handpieces

		4 (Drill depth, mm)	6 (Drill depth, mm)	
Color bars reflect family color shown on the adapter	8-flute	 <p>PowerPick Attachment Designed to make quick work of microfracture procedures by using the shaver handpiece to power the Ø1.5 mm drill tip; 30° and 45° angled shafts facilitate use in a variety of other applications, including ACL femoral tunnel location and crimson duvet techniques in the shoulder</p>	AR-8150PP-30 (30°) AR-8150PP-45 (45°)	AR-8150PX-45 (45°)
	8-flute	 <p>PowerRasp Attachment Flat, reciprocating device that creates a flat, smooth bony surface and helps prevent tissue wrapping; best suited for acromioplasties, distal clavicle resections, and osteophyte resections</p>	AR-8350PR	AR-8400PR AR-8550PR





NanoResection™ Small-Hub Shaver Handpiece Options (7 cm)

		1 (Ø mm)	2 (Ø mm)	3 (Ø mm)	3.5 (Ø mm)	
Color bars reflect family color shown on the adapter		PowerRasp™ Resector, Small Hub Flat, reciprocating device to create a flat, smooth bony surface and help prevent tissue wrapping; ideally suited for osteophyte resection or bony site preparation		AR-9350PR		
		PowerPick™ Drill, Small Hub Excellent for microdrilling subchondral bone defects		AR-9100PP-00 (0°) AR-9100PP-45 (45°)		
		Torpedo™ Shaver Blade Tapered tip and scissor-like cutting action facilitate rapid, aggressive soft-tissue resection		AR-9350TD		
		Dissector, Small Hub Sharp, teardrop-shaped outer window; toothed inner cutting window aggressively attacks tissue, grabbing it and drawing it into the cutting area		AR-9300DS	AR-9350DS	
		Sabre, Small Hub Sharp, teardrop-shaped outer window; smooth inner cutting window leaves crisp, clean edges for a more anatomic appearance following resection		AR-9200SR	AR-9300SR	AR-9350SR
	10-flute		Oval Burr, Small Hub Designed for rapid, precise bone resection for small joint applications		AR-9300OBT	AR-9350OBT
	10-flute		Round Burr, Small Hub A precise and aggressive burr designed for osteophyte removal in small joint applications		AR-9300RBT	AR-9350RBT

NanoResection™ Instruments (11 cm length)

Fits Small-Hub Shaver Handpiece














		2.8 (Ø mm)	
10-flute		Nano Sabre Shaver Unique 2.8 mm × 11 cm size provides length to reach all necessary pathology while allowing for "sutureless incisions" and improved navigation around tight joint spaces; sharp, teardrop-shaped outer window; smooth inner cutting window leaves crisp, clean edges for a more anatomic appearance	AR-9280NSR
		Nano Oval Burr Nano-sized oval burr designed for precise bone removal while supporting "sutureless incisions"; 10 extrasharp flutes for improved efficiency; distal cutting capabilities for tight joint spaces	AR-9280NOT

ShaverDrill™ Instruments

Fits Large-Hub Shaver Handpieces

ShaverDrill instruments quickly attach directly to our large-hub shaver handpieces and are used to predrill pilot holes for a variety of anchors. By eliminating the need to open a battery-powered drill tray, ShaverDrill instruments

may help to extend battery life, reduce intraoperative setup time, and lower cleaning and processing costs between cases.

Color bars reflect family color shown on the adapter		ShaverDrill Instrument Kit for 1.6 mm FiberTak® anchor	AR-3600NDSR-1
		Knotless FiberTak drill, 1.8 mm	AR-3610NSD-2
		ShaverDrill instrument for 2.4 mm PushLock anchors	AR-2922DSR-24
		ShaverDrill instrument for 2.4 mm PushLock anchors, hard bone	AR-2922DTSR-24
		ShaverDrill instrument for 2.6 mm Knotless FiberTak soft anchor	AR-3657SD
		ShaverDrill instrument for 2.9 mm PushLock® anchors	AR-2923DSSK
		ShaverDrill instrument for short 2.9 mm PushLock anchors	AR-2923DSR
		ShaverDrill instrument for short 2.9 mm PushLock anchors, hard bone	AR-2923DTSR
		ShaverDrill Kit for short 2.9 mm PushLock anchors, hard bone	AR-2923DTSK
		ShaverDrill instrument for 3 mm SutureTak® anchors	AR-1250LTSR
		ShaverDrill Kit for 3 mm SutureTak anchors	AR-1934DSSK
		ShaverDrill instrument for 3 mm SutureTak anchors, hard bone, and Knotless SutureTak anchors	AR-1938DSR
		ShaverDrill Kit for 3 mm SutureTak anchors, hard bone, and 3 mm Knotless SutureTak anchors	AR-1938DSSK

Fluid and Insufflation Management

68	DualWave™ Arthroscopy Pump
69	Continuous Wave™ 4 Arthroscopy Pump and Accessories
70	Main Pump Inflow Tubing
70	ReDeuce™ 2-Piece Pump Inflow Tubing System
70	Extension Inflow Tubing
70	DualWave™ Pump Tubing
71	Inflow/Outflow Tubing Set
71	Gravity Tubing
72	Waste Collection Bag
72	Y-Tubing Adapter
72	Autoclavable Remotes
72	Footswitch
72	Cannulas and Parts
73	DualWave™ Cart
73	Continuous Wave™ 4 Cart
74	Synergy ^{Insufflation} Console
75	The Suctioner® Surgical Suction Mat
75	The EcoSuctioner® Surgical Suction Mat
76	The SuctionBlu™ Mat
76	The Surgical Mat™ and The Mini Mat™ Surgical Mats
77	EcoDri-Safe™ Absorbent Rolls
77	Dri-Safe® Absorbent Pads
78	Fluid Diamond® Floor Aspirator

DualWave™ Arthroscopy Pump



The DualWave integrated inflow and outflow fluid management system may also be used as a simple inflow-only arthroscopy pump. The DualWave arthroscopy pump can solve any fluid management need, whether it is controlling distention during complex arthroscopic procedures or helping a facility increase cost-containment capabilities.

Clear, consistent, turbulence-free visualization is achieved with the DualWave arthroscopy pump using a proprietary design that independently controls the shaver and cannula outflow to maintain optimal joint distention and clarity.

The touch panel video display gives real-time pressure and flow readings and allows the user to track the total fluid used during a case and the total operating time at the end of a procedure. User-programmable defaults encompass pressure presets (by joint space), outflow rate, and shaver suction levels, giving the surgeon a consistent, repeatable operating experience.

Flexibility is key with touch panel access and the choice of a foot pedal or an autoclavable remote control that allows the user to control the DualWave arthroscopy pump from the sterile field. The Rinse and Lavage modes allow the user to control all aspects of fluid management at the touch of a button. Lavage mode provides a timed cycle of increased pressure to help maintain proper visualization, and the Rinse mode increases the fluid passing through the joint space without affecting the pressure. When operating in Inflow Only and Inflow and Outflow modes, the shaver boost function automatically compensates for the loss of distention that can occur when shaver suction is introduced into the joint space by increasing the pressure while the shaver is active.

To help control costs, the DualWave arthroscopy pump can use the ReDeuce tubing option to significantly lower the tubing and fluid cost per patient, and lead to reduced room turnover times.

DualWave arthroscopy pump

AR-6480

Continuous Wave™ 4 Arthroscopy Pump and Accessories



Arthrex offers a clear vision of the future with the Continuous Wave 4 arthroscopy pump. The Continuous Wave 4 pump, now with a Synergy heads-up display and direct communication with the Synergy^{Resection}™ system, is a fully integrated inflow-only arthroscopy pump.

Clear, consistent, turbulent-free visualization is achieved through Arthrex's latest proprietary design, which maintains optimal joint distention and clarity.

The touch-panel video display gives real-time pressure and flow readings, and displays the total fluid used during a case, as well as the total run time at the end of a procedure. Flexibility is key with touch-panel access and the choice of a foot pedal or autoclavable remote-control device that allows the user to control the Continuous Wave 4 arthroscopy pump from the sterile field.


User-selected defaults encompass pressure settings (by joint space), flow rate, and shaver boost levels. The lavage mode allows the user to control pressure at the touch of a button. The addition of shaver detection provides an optional boost of pressure when the shaver is activated.

Continuous Wave 4 arthroscopy pump

AR-6485

Inflow Pump Tubing and Outflow Pump Tubing

Main Pump Inflow Tubing



This set is for those who prefer a complete tubing set and want to replace the fluid bags after each case. Each tubing set comes with an independently packaged press-fit Luer lock to allow for tubing length adjustment if used in conjunction with the extension tubing. If used alone, the set must be discarded after each case.


Main pump tubing for AR-6480 or AR-6485, 13.5 ft long, sterile, qty. 10 AR-6410

Extension Tubing



This unique tubing provides an economical option to reduce tubing costs by allowing the main pump tubing to remain in use for the entire surgical day. The exclusive backflow valve prevents contaminated fluid from reaching the main pump tubing. The main pump tubing and fluid bags no longer need to be discarded after each surgical case.

Extension tubing, 8 ft long, sterile, qty. 20 AR-6220



DualWave outflow tubing allows the DualWave arthroscopy pump to regulate both inflow and outflow in the joint space. By controlling both aspects, distention is maintained at the highest level, resulting in maximum visualization.

DualWave outflow tubing for AR-6480, qty. 10 AR-6430

ReDeuce™ 2-Piece Pump Inflow Tubing System




The 2-piece pump tubing maintains sterility due to the design of the patient-side tubing and the touch-proof connector. Leaving the pump-side tubing in place the entire day saves on fluid because the remaining fluid does not need to be discarded after terminating each surgical case.

ReDeuce 2-piece pump tubing w/ connector for AR-6480 or AR-6485, 8 ft long, sterile, qty. 10 AR-6411



The 2-piece patient tubing set contains a backflow restrictor, saving tubing costs and setup time by only requiring patient tubing to be discarded and passed off from the sterile field after and before each case.

ReDeuce 2-piece patient tubing, sterile, 8 ft long, qty. 20 AR-6421



Tube set, inflow, arthroscopy AR-6413

Inflow/Outflow Combination Tube Sets and Gravity Tubing

Inflow/Outflow Combination Tube Sets



This combination set comes with both AR-6410 and AR-6430, simplifying preparation by reducing the number of items needed in the OR to set up the DualWave™ pump. Fewer tubing packages reduce ordering and simplify inventory tracking, providing a direct price reduction.

DualWave outflow tubing set w/ main pump tubing for AR-6480, qty. 8	AR-6436
---	---------



This combination set comes with both AR-6421 and AR-6430, simplifying preparation by reducing the number of items needed in the OR to set up the DualWave pump. Fewer tubing packages reduce ordering and simplify inventory tracking, providing a direct price reduction.

DualWave 2-piece pump tubing w/ connector for AR-6480, qty. 10	AR-6435
---	---------

Gravity Tubing



2-Spike



4-Spike

Arthrex Gravity Tubing is compatible with all arthroscopic sheaths and cannulas and able to maintain proper distention even in complicated surgical cases that demand high flow.

2-spike tubing	AR-6412
4-spike tubing	AR-6414

Pump Accessories

Waste Collection Bag



The waste collection bag can be hung from the AR-6481 pump cart and can hold up to 10 liters of waste fluid. At the end of a surgical case, simply seal the bag and dispose.

Waste collection bag, qty. 10 AR-6431

Autoclavable Remote



The autoclavable DualWave™ remote control provides finger-touch access to pressure adjustment, Lavage and Rinse modes, shaver suction levels, and the Run/Stop button.

DualWave autoclavable remote AR-6495

Cannulas and Parts



Flexible cannula,
2.5 mm × 103 mm with
obturator

AR-3420.2.5

Y-Tubing Adapter



The optional Y-tubing adapter allows 4 fluid bags to be connected to either the main pump tubing or the 2-piece pump tubing, resulting in fewer fluid bag changes.


Y-tubing adapter, double spike, 2
per set, sterile, qty. 10 AR-6215

Footswitch



The DualWave footswitch allows the user to control both pressure and flow levels with the lavage and rinse pedals.

DualWave footswitch AR-6483



Flexible cannula,
4.5 mm × 75 mm **(a)**

AR-6525

Valve for cannula **(b)**

AR-6501

Trocar for cannula


AR-6502

Obturator for cannula **(c)**

AR-6503

Pump Accessories: Cart

DualWave™ Cart



The DualWave arthroscopy pump cart provides easy installation of the DualWave or the Continuous Wave 4 pump at the proper height for surgery. The mobile cart can be configured for storage of fluid bags, shaver blades, and/or foot pedals for easy rotation between operating rooms.

DualWave cart

AR-6481

Continuous Wave™ 4 Cart



The Continuous Wave 4 pump cart provides easy installation at the proper height for surgery. The mobile cart features 4 saline bag hooks, 1 pump stand, and 1 storage basket.

Continuous Wave 4 cart

AR-6485PC

Synergy^{Insufflation}™ System

Synergy^{Insufflation} Console



(a)

The technologically sophisticated yet simple Synergy^{Insufflation} device is used for diagnostic and therapeutic endoscopic procedures to distend the peritoneal cavity by filling it with gas. The device has a volumetric flow rate of up to 50 L to quickly fill the abdomen. The intuitive touchscreen interface allows for quick and easy adjustments. Pressure and volumetric flow rate can be adjusted between High-Flow and Bariatric modes according to preference. A built-in heating element gives surgeons the ability to heat the gas before it is delivered into the patient. Using our Synergy video management system, the heads-up display feature allows for easy visibility of system settings on the monitor. The console features a venting system and a contamination alarm.

Synergy ^{Insufflation} console (English, German, French, Spanish) (a)	AR-3290-0004
Insufflation tube set with ISO connection, disposable, sterile, 20 pcs./box	Z0293
Tube set ISO connector, reusable 20×, autoclavable	AR-3290-1002
Disposable heated tube set (pack of 10)	AR-3290-1007
Heated tube set, reusable 100×, autoclavable	AR-3290-1008
Dual-yoke CO ₂ switching assembly	002433-00
ISO barb adapter	Z0164-01
High-pressure hose US bottle/PIN index, 1.5 m	Z0499-01
High-pressure hose US bottle/DIN, 1.5 m	Z5044-01
High-pressure hose US bottle/ISO, 1.5 m	Z5045-01
Sinter filter for gas connection	Z5072-01
Gas connection, US	Z5073-01
Gas connection, NIST	Z5074-01
Gas connection, DISS	Z5075-01
Gas connection, universal, angled	Z5077-01
Gas connection, universal	Z5078-01
US/DIN hose for house gas, 3 m	Z5100-01
US/DIN hose for house gas, 5 m	Z5101-01
US/NORME French hose for house gas, 3 m	Z5102-01
US/NORME French hose for house gas, 5 m	Z5103-01

Disposable Antifatigue Mats

The Suctioner® Surgical Suction Mat



This disposable surgical mat combines the absorbency of the poly-backed Dri-Safe® pad and the nonslip, anti-fatigue memory foam construction of The Surgical Mat™ floor covering with powerful, quiet, and continuous suction capability. These features are designed to ensure a drier, safer, and cleaner operating room. The Suctioner mat has graduated slots to move fluid toward lateral suction tubes.

This mat easily attaches to any OR vacuum source in seconds and efficiently manages unwanted fluids, keeping them off the OR floor. The mat contains antimicrobial properties designed to reduce the risk of cross-contamination. This latex-free surgical mat is ideal for all surgical procedures.

The Suctioner surgical mat, 22 in × 35 in (comfort mat), USA-X10
42 in × 56 in (total), blue, qty. 8

The EcoSuctioner® Surgical Suction Mat



- › All-in-one ergonomic suction mat that allows 2 people to stand comfortably during a procedure
- › Designed with directional grooves for improved quiet, continual suction; attaches to any suction source
- › Large, rugged absorbent mat included on underside creates more nonslip surface in OR
- › Lightweight suction mat, weighs 0.79 kg, allowing for quick and easy setup and removal
- › Reduces laundry on the floor

The EcoSuctioner surgical mat, 28 in × 40 in (comfort mat), USA-XG8
42 in × 56 in (total), green, qty. 8

The Suctioner, Dri-Safe, The Surgical Mat, and The EcoSuctioner are trademarks of Flagship Surgical, LLC.

The SuctionBlu™ Mat



The SuctionBlu disposable surgical mat combines the absorbency of the poly-backed Dri-Safe® pad and the nonslip, anti-fatigue memory foam construction of The EcoSuctioner® mat with powerful, quiet, continuous suction capability. These features are designed to create a drier, safer, and cleaner operating room.


The graduated slots on The SuctionBlu mat use gravity to direct fluids toward the suction tube. The mat easily attaches to any OR vacuum source in seconds and efficiently manages unwanted fluids, keeping them off the floor.

The SuctionBlu mat is a smaller, more economical version of The EcoSuctioner mat but larger than The Suctioner® mat.

The SuctionBlu surgical mat, 42 in × 56 in, blue

USA-XBL

The Surgical Mat™ and The Mini Mat® Surgical Mats



The Surgical Mat and The Mini Mat surgical mats were designed to address three OR issues: comfort, safety, and fluid buildup. Both mats are constructed of a resilient memory foam material that provides relief and support when standing on hard concrete floors and steel step stools. The mats feature a nonslip molded foam design and custom tapered fluid channels to keep the surgical team's feet dry and reduce the risk of tracking body fluids across the OR. The Surgical Mat and The Mini Mat surgical mats are latex-free and ideal for all surgical procedures.

The Surgical Mat pad, 20 in × 39 in, blue, qty. 10 **(a)**

USA-010

The Mini Mat surgical mat, 13 in × 17 in, blue, qty. 10 **(b)**

USA-A10

The SuctionBlu, Dri-Safe, The EcoSuctioner, The Suctioner, The Surgical Mat, and The Mini Mat are trademarks of Flagship Surgical, LLC.

Fluid Containment Pads

EcoDri-Safe™ Absorbent Rolls



Containing fluids and streamlining clean-up procedures has never been easier than with the EcoDri-Safe absorbent rolls. The 28 in × 100 ft rolls are perforated every 20 in to create 60 easy-to-tear-off 20 in × 28 in floor mats on every roll.

EcoDri-Safe absorbent rolls are ideal wherever fluid containment can be an issue:

- › Operating rooms
- › Scrub sinks
- › Laboratory stations
- › Central processing
- › Emergency rooms
- › Procedure rooms
- › Labor and delivery wards

EcoDri-Safe non-poly-backed roll, 100 ft × 28 in, 20 in × 28 in, 60 sheets per roll	USA-ECO100
EcoDri-Safe poly-backed roll, 100 ft × 28 in, 20 in × 28 in, 60 sheets per roll	USA-ECO100PLY
EcoDri-Safe dispensing cart	USA-ECODISP

Dri-Safe® Absorbent Pads




- › Superabsorbent, economical, disposable floor pads that eliminate expensive laundry costs
- › Use wherever fluids build up, including OR floors, in front of scrub sinks, lab station floors, etc
- › Dri-Safe pads are designed to help prevent slips and falls and the transmission of harmful contaminants to health care providers
- › Impermeable poly-backing option promotes safety by keeping the pads firmly in place and locking in fluids

Dri-Safe absorbent pad, 28 in × 40 in, blue, qty. 20	USA-B28X40
Dri-Safe absorbent pad, 28 in × 40 in, blue, poly-backing, qty. 20	USA-B28X40PLY
Dri-Safe absorbent pad, 28 in × 56 in, blue, qty. 15	USA-B28X56
Dri-Safe absorbent pad, 28 in × 56 in, blue, poly-backing, qty. 15	USA-B28X56PLY

EcoDri-Safe and Dri-Safe are trademarks of Flagship Surgical, LLC.

Floor Aspirator

Fluid Diamond® Floor Aspirator



- › Fast, effective fluid collection
- › Powerful suction action designed to keep floors dry and reduce the potential for slipping or falling
- › Eliminates the need for slippery, costly towels and blankets on the OR floor to absorb fluids
- › Designed with quieting pad and anti-floor-sticking technology
- › Attaches to any surgical suction source
- › Includes easily attachable suction tubing

Fluid Diamond floor aspirator, qty. 30

USA-FDX

Fluid Diamond is a trademark of Flagship Surgical, LLC.



Arthroscopic RF System

82	Synergy ^{RF™} Console
83	Synergy ^{RF™} Footswitch
83	Apollo ^{RF®} i90 Probe
84	Apollo ^{RF®} MP90 Probe and Apollo ^{RF} XL90 Probe
85	Apollo ^{RF®} Hook Probe
85	Apollo ^{RF} MP50 Probe
86	Apollo ^{RF®} H50 Probe
86	Apollo ^{RF} SJ50 Probe

Synergy^{RF}™ Bipolar Ablation System

The Synergy^{RF} console is an advanced electro-surgical generator designed to meet the needs of arthroscopic surgical procedures. The console's advanced output control algorithm uses feedback systems that are capable of adjusting the output as the tissue characteristics change. This will result in a more consistent clinical effect at the surgical site. The Synergy^{RF} console provides an interface to the Synergy Video Management System for a heads-up display feature so settings are displayed on the main surgical monitor. The Synergy 4K Imaging System allows the user to apply surgeon- and procedure-specific presets for ablation and coagulation. These presets can be applied for each Apollo^{RF} probe in the system.

With the ablation probes, Arthrex offers a complete bipolar RF system for the resection, ablation, and coagulation of soft tissue and hemostasis of blood vessels in arthroscopic and orthopedic procedures.

The system features intelligent device recognition (IDR) for the probes currently in use. The Metal Proximity Detection function suspends power when the probe is close to the metal portion of the scope or metal instruments. The ergonomic handle of the Apollo^{RF} probe features three buttons to operate the ablation and coagulation functions and to adjust the Ablate power setting.

Synergy^{RF} Console



- › Integrates with the Synergy Imaging System or operates as a standalone unit
- › Probe Recognition function identifies probe currently in use
- › Metal Proximity Detection function suspends power when in close proximity to the metal portion of the scope or metal instruments
- › User-friendly LCD touchscreen
- › Easy-to-adjust ablation and coagulation settings
- › Uses controlled plasma-based bipolar RF technology with Apollo^{RF} probes

Synergy^{RF} console

AR-9800

Synergy^{RF} Footswitch



- › Separate Ablate and Coag pedals
- › Central ablation-level-adjustment button
- › Integrated carrying handle

Synergy^{RF} foot pedal

AR-9800-F

Apollo^{RF} i90 Probe



- › Torpedo-shaped electrode* designed for easy access into the joint and accurate removal of tissue
- › Increased flow rate/improved flow path* designed to efficiently remove bubbles for better visual clarity, reduced tissue clogging, and faster removal of fluid from the joint¹
- › Edge Control feature on the electrode provides greater precision and gentle ablation distally and on the sides of the electrode
- › Enhanced ablation performance for efficient tissue removal
- › Shorter working length of 140 mm

Apollo^{RF} i90 probe


AR-9831

Reference

1. Arthrex, Inc. Data on file (AR-9831 clog performance design verification). Naples, FL; 2021.

*Patent pending

Apollo^{RF} MP90 Probe and Apollo^{RF} XL90 Probe



- > 90° tip configured to access and ablate tissue rapidly
- > Durable one-piece active electrode
- > Easy-to-use, ergonomic hand controls
- > Multiple ports for efficient bubble removal

- > Apollo^{RF} XL90 probe has a 25% longer electrode than the Apollo^{RF} MP90 probe to remove more tissue per pass

Apollo ^{RF} MP90 probe	AR-9811
Apollo ^{RF} XL90 probe	AR-9821

Apollo^{RF} Hook Probe



- › Generous hook size for easy tissue capture and rapid tissue resection and coagulation
- › Durable one-piece active electrode
- › Easy-to-use, ergonomic hand controls

- › Configured to access and reach tissue in various arthroscopic procedures

Apollo^{RF} hook probe

AR-9825

Apollo^{RF} MP50 Probe



- › Smooth radius at tip designed to access confined anatomy
- › 360° edge control for precise ablation
- › Dual aspiration slots for effective suction

- › Lower ablation default setting of 4 for use around sensitive tissue
- › For use in knee and shoulder procedures

Apollo^{RF} MP50 probe

AR-9815

Apollo^{RF} H50 Probe



- › 185 mm working length provides an extended reach into the hip joint
- › Reinforced shaft and handle junction maintain device integrity during hip arthroscopy
- › 360° edge control for precise ablation and debridement
- › Flat electrode face for smooth contouring of tissues
- › Optimal tip geometry for enhanced access of hip anatomy
- › Dual aspiration slots for effective suction
- › Upgraded insulation material to increase durability during portal entry

Apollo^{RF} H50 probe

AR-9835

Apollo^{RF} SJ50 Probe



- › Reduced probe dimension of 110 mm working length and 3.3 mm diameter for ease of access in confined anatomy
- › 360° edge control for precise ablation
- › Anatomic 50° tip configured to reach and ablate tissue efficiently
- › Low default ablation setting of 4 for reduced energy delivery
- › Optimal for foot, ankle, elbow, knee, and Nano arthroscopy procedures
- › Optimized aspiration to remove tissue particulate and maintain optimal fluid temperatures

Apollo^{RF} SJ50 probe

AR-9845



Synergy Service and Support

Service and Support



Arthrex Armour Technology Protection

Service You Can Rely On

We understand that patients rely on you every day—and that you rely on your equipment. You have purchased the most innovative world-class imaging devices available, and you always expect them to be ready.

As your strategic partner, you can rely on Arthrex to provide quality products and service to continue Helping Surgeons Treat Their Patients Better®.

Arthrex Armour offers customizable service plan levels to meet all your product and service needs.

Arthrex Armour Core Protection

Arthrex Armour Core Protection provides our baseline service level support billed per-incident and is available on all serviceable Arthrex products. Prepaid per-incidents are also available for all scope product lines, Synergy Vision™ consoles and camera heads, Synergy^{HD}™ consoles and camera heads, and Synergy^{UHD4}™ consoles and camera heads. With prepaid per-incidents, you can better forecast your operational budget through a bulk purchase.

This service level includes all the following benefits:

- › Advanced replacements
- › Damage protection
- › Free replacement shipping
- › 24-hour remote technical support

Arthrex Armour Premium Protection

This level of Arthrex Armour provides total protection against product issues and incidental damage.

Multiple payment options make an Arthrex Armour Premium Protection plan a budget-friendly solution that guarantees peak equipment performance and helps ensure maximum uptime.

This service level includes all the following benefits:

- › Advanced replacements
- › Damage protection
- › Free replacement shipping
- › Free return shipping
- › 24-hour remote technical support and software updates
- › Onsite technical support

Arthrex Armour Titanium Protection

To resolve any possible technical challenges, Arthrex Armour Titanium Protection offers our most comprehensive service. Along with expert 24-hour remote technical support provided by our Technical Assistance Center, this tier features a trained Arthrex agency representative residing on location to provide immediate technical support. Layer this protection plan with Arthrex Armour Premium Protection to ensure the best service experience possible.

This service level includes all the following benefits:

- › 24-hour remote technical support
- › Onsite Arthrex agency representative for immediate support
- › Onsite software updates

Why Arthrex?

Our commitment to using high-quality, certified components means we can provide the highest standard of service designed to keep your equipment in excellent working condition over time. For maximum uptime and the highest-quality performance, we are the logical choice for service replacements and after-sale care.

Contact Information

Please contact us at ServiceAgreements@arthrex.com with any questions or requests regarding our service offerings.

Warranty Policy

For all warranty information, including disclaimers, exclusions, terms, conditions, and related provisions, refer to the Arthrex U.S. Product Warranty section of the Arthrex, Inc., website, found at www.arthrex.com/corporate/arthrex-us-product-warranty.

This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.



Arthrex manufacturer, authorized representative, and importer information (Arthrex eIFUs)



US patent information

