



OrthoLine™ System Fracture Plate

Size Chart for Distal Radial Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Radial Plate, Ti, 1.6 mm		VAR-3116DR-L	39.0	2.6	1.4	7.3	5.0
		VAR-3116DR-R					
Distal Radial Plate, Ti, Broad, 1.6 mm		VAR-3116BDR-L	42.3	3.1	1.8	7.9	5.4
		VAR-3116BDR-R					
Distal Radial Plate, Ti, 2.0 mm		VAR-3120DR-L	46.9	3.7	1.9	8.8	5.7
		VAR-3120DR-R					
Distal Radial Plate, Ti, 2.4 mm		VAR-3124DR-L	57.2	4.7	2.5	10.0	7.1
		VAR-3124DR-R					



OrthoLine™ System Fracture Plate

Size Chart for Distal Radial Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Radial Plate, Ti, 3.0 mm		VAR-3130DR-L	64.1	5.4	3.2	11.6	8.8
		VAR-3130DR-R					
Distal Radial Plate, SS, 3.5 mm		VAR-3035DR-L	80.1	6.3	4.1	13.6	11.4
		VAR-3035DR-R					

