OrthoLine™ Fracture Management System

Weight Reference Chart: Straight Plates, T-Plates, and Cuttable Plates

	₩	1	H	×	H		
	1-4 kg	4-8 kg	8-12 kg	13-19 kg	20-30 kg	31-45 kg	45-55+ kg
Radius	1.6 mm, 1.6 mm B	2.0 mm	2.0 mm, 2.4 mm	3.0 mm	3.0 mm, 3.5 mm	3.5 mm	3.5 mm
Humerus	1.6 mm, 1.6 mm B	2.0 mm, 2.4 mm	2.4 mm, 3.0 mm	3.0 mm, 3.5 mm	3.5 mm	3.5 mnm, 3.5 mm B	3.5 mm B
Femur	1.6 mm, 1.6 mm B	2.0 mm, 2.4 mm	2.4 mm, 3.0 mm	3.0 mm, 3.5 mm	3.5 mm	3.5 mm, 3.5 mm B	3.5 mm B
Tibia	1.6 mm, 1.6 mm B	2.0 mm, 2.4 mm	2.4 mm	3.0 mm, 3.5 mm	3.5 mm	3.5 mm, 3.5 mm B	3.5 mm B
T-Plate	1.6 mm, 1.6 mm B	2.0 mm	2.0 mm, 2.4 mm	3.0 mm	3.5 mm	3.5 mm	3.5 mm

All weight ranges are suggestions based on the ultimate load of the product and/or our veterinary consultant's recommendation. It is up to the veterinarian to determine the proper product, application, technique, and weight range for its use.

