

OrthoLine™ Distal Humeral Fracture Management System

Weight Reference Chart



	1 kg-6 kg	6 kg-8 kg	8 kg-13 kg	13 kg-20 kg	20 kg-35 kg	
Humerus	2.0 mm	2.0 mm	2.4 mm	3.0 mm	3.5 mm	
Screw Range	1.6 mm/2.0 mm	2.0 mm	2.0 mm/2.4 mm	2.4 mm/3.0 mm	2.7 mm/3.5 mm	
Plate						
						

All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultants' recommendations. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.