










OrthoLine™ System Fracture Plate







Size Chart for Proximal Femur Plate

	Product Image	Item Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Proximal Femur Plate, Ti, Broad, 1.6 mm		VAR-3116BPF-L	64.2	7.0	1.9	7.0	5.0
		VAR-3116BPF-R					
Proximal Femur Plate, Ti, 2.0 mm		VAR-3120PF-L	78.6	7.9	2.2	8.1	6.1
		VAR-3120PF-R					
Proximal Femur Plate, Ti, 2.4 mm		VAR-3124PF-L	96.5	9.7	2.9	9.5	7.6
		VAR-3124PF-R					
Proximal Femur Plate, Ti, 3.0 mm		VAR-3130PF-L	154.3	11.5	3.6	12.6	9.4
		VAR-3130PF-R					



OrthoLine™ System Fracture Plate

Size Chart for Proximal Femur Plate

	Product Image	Item Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Proximal Femur Plate, SS, Short, 3.5 mm		VAR-3035PFS-L	155.1	13.0	4.3	15.2	11.9
		VAR-3035PFS-R					
Proximal Femur Plate, SS, 3.5 mm		VAR-3035PF-L	188.1	13.0	4.3	15.2	11.9
		VAR-3035PF-R					
Proximal Femur Plate, SS, Broad, 3.5 mm		VAR-3035BPF-L	201.7	14.3	5.5	17.7	14.0
		VAR-3035BPF-R					