




OrthoLine™ Fracture Management System

Weight Reference Chart – Proximal Femoral Plate



	1 kg-4 kg	4 kg-6 kg	7 kg-12 kg	13 kg-19 kg	20 kg-32 kg	31 kg-46 kg	45 kg-55 kg
Plate Size	1.6B mm	1.6B mm/2.0 mm	2.4 mm	3.0 mm	3.5 mm	3.5 mm/3.5B mm	3.5B mm
Plate							

All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultants' recommendations. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.