

OrthoLine™ Fracture Management System—Iliac Fracture Plate

Weight Reference Chart



	1-4 kg	3-9 kg	9-14 kg	14-25 kg	23-50 kg
Plate Size	1.6 mm	1.6 mm/2.0 mm	2.4 mm	3.0 mm	3.5 mm
Plate					

All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultants' recommendations. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.