

# | Arthrex TPLO Saw System



**Arthrex®**  
Vet Systems

# Handpiece and Accessories

Available as a complete kit or as individual components

## Key Features and Benefits

- › Linear form factor dedicated TPLO handpiece
- › Nickel-metal hydride batteries (2)
- › Charging adapter
- › Aseptic transfer sleeves (2)
- › Allen wrench
- › System transport case
- › Sterilization tray



# TPLO Blade Options

Available individually or in cost-savings sets



Fenestrated	8 mm	10 mm	12 mm	15 mm	18 mm	21 mm	24 mm	27 mm	30 mm	33 mm
Nonfenestrated		10 mm	12 mm	15 mm	18 mm	21 mm	24 mm	27 mm	30 mm	33 mm

# Arthrex Armour

## Technology Protection

- › Arthrex Armour offers service plans to meet all your product and service needs.
- › Contact your local Arthrex Vet Systems representative for more information.



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience, and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.



Arthrex manufacturer, authorized representative, and importer information (Arthrex eIFUs)



US patent information