


1.6 mm OrthoLine™ System


VAR-3116S-01



Pic.	Item Number	Qty.	Description
1	VAR-4020-01	1	Screwdriver, T6 hexalobe, self-retaining
2	VAR-4020-02	1	Driver, T6 hexalobe, self-retaining
3	VAR-2024DD	1	Depth device, 1.6 mm/2.0 mm/2.4 mm
4	VAR-4016TDG	1	Drill guide, 1.1 mm
5	VAR-4016VDG	1	Drill guide, variable, 1.6 mm
6	VAR-4016DG	1	Drill/depth guide, locking, 1.6 mm
7	VAR-4016D	2	Drill bit, solid, AO, 1.1 mm
8	VAR-4016SD	2	Drill bit short, solid, AO, 1.1 mm
9	VAR-4000GC-01	1	General sterilization case insert
10	VAR-8941F	1	Screw-holding forceps
11	VAR-3016SC-01	1	1.6 mm screw caddy, VAL
12	VAR-4000GC	1	General sterilization case

1.6 mm OrthoLine™ System

VAR-3116S-01



Pic.	Item Number	Qty.	Description
13	VAR-4000-07	2	Bending iron, small
14	VAR-4000BPC	1	Bending plug caddy
15	VAR-4020-04	5	Bending plug, cannulated, 1.6 mm/2.0 mm

*BB-Taks are not included in set

1.6 mm Screw Set

VAR-3116SCS-01



Pic.	Item Number	Qty.	Description
1	VAR-8916-06-20	4 ea	Low-profile cortical screw, ti, 6-20 mm
2	VAR-8916V-06-20	4 ea	Low-profile variable screw, ti, 6-20 mm

*Screw caddy and screw-holding forceps are not included in screw set

1.6 mm OrthoLine™ System

VAR-3116SCS-01

Date of Surgery:

Surgeon Name:

Account:

Representative Contact Information:

Set Total:

Miscellaneous Notes:

To obtain a copy of a product's directions for use (DFU), which includes sterilization instructions, please visit our website at www.arthrex.com/edfu or contact your local representative.

arthrexvetsystems.com

© 2026-07 Arthrex, Inc. All rights reserved. vDOC1-002271-en-US_A

