

Knotless SwiveLock[®] Anchor System

Instrumentation Charts



Arthrex[®] 
Vet Systems

3.5 mm Knotless SwiveLock Kit

VAR-2325PSLK

	Product Description	Item Number	Product Image
Small Breed Weight Range 16-40 lb	Drill Guide/Tissue Protector for 3.5 mm SwiveLock Anchor	VAR-1678-01	
	Drill for 3.5 mm SwiveLock Anchor	VAR-1678-05	
	Tap for 3.5 mm SwiveLock Anchor	VAR-1678-03	
	Cannulated Drill Bit, 2.0 mm (AO quick connect)	VAR-8933-20C	
	Guidewire, 0.045 in (1.14 mm)	VAR-8933K	
	Suture Passing Wire, Nitinol, 8 in (20.3 cm)	VAR-1255-08	
	Small Joint FiberWire® Scissors	VAR-11797	
	CrCL Femoral Socket Aiming Device	VAR-8920T	

Weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultant's recommendation. It is up to the veterinarian to determine the proper product, application, technique, and weight range to use the selected product.



4.75 mm Knotless SwiveLock Kit

VAR-2324PSLK

	Product Description	Item Number	Product Image
Medium Breed Weight Range 30-65 lb	Spade Tip Drill for 4.75 mm and 5.5 mm SwiveLock Anchor	VAR-1927D	
	Tap for 4.75 mm SwiveLock Anchor	VAR-2324PTB	
	Cannulated Drill Bit, 2.5 mm (AO quick connect)	VAR-8737-09	
	Guidewire, 0.045 in (1.14 mm)	VAR-8933K	
	Suture Passing Wire, Nitinol, 8 in (20.3 cm)	VAR-1255-08	
	Small Joint FiberWire® Scissors	VAR-11797	
	CrCL Femoral Socket Aiming Device	VAR-8920T	

Weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultant's recommendation. It is up to the veterinarian to determine the proper product, application, technique, and weight range to use the selected product.



5.5 mm Knotless SwiveLock Kit

VAR-2323PSLK

	Product Description	Item Number	Product Image
Large Breed Weight Range 40-80 lb (80-120 lb requires 2 pieces of FiberTape® sutures)	Spade Tip Drill for 4.75 mm and 5.5 mm SwiveLock Anchor	VAR-1927D	
	Tap for 5.5 mm SwiveLock Anchor	VAR-1927CTB	
	Cannulated Drill Bit, 2.5 mm (AO quick connect)	VAR-8737-09	
	Guidewire, 0.045 in (1.14 mm)	VAR-8933K	
	Suture Passing Wire, Nitinol, 8 in (20.3 cm)	VAR-1255-08	
	Small Joint FiberWire® Scissors	VAR-11797	
	CrCL Femoral Socket Aiming Device	VAR-8920T	

Weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultant's recommendation. It is up to the veterinarian to determine the proper product, application, technique, and weight range to use the selected product.

