

Drill System Specifications

Reference Chart

Product Image	Drill	Handpiece Weight	Drill Speed	Drill Torque	Reamer Speed	Reamer Torque	Cutting Speed
	DrillSaw Mini 300™ system	550 g (1.2 lb)	1300 rpm	24 N·cm (2 in·lb)	500 rpm	1.15 N·m (10 in·lb)	15,000 rpm
	DrillSaw Sports 400™ system	630 g (1.3 lb)	1300 rpm	3.2 N·m (28 in·lb)	330 rpm	9 N·m (79 in·lb)	22,000 rpm