




















Arthrex TPLO Saw Blade

Weight/Radius/Plate Size Reference Chart

										
	1.6 mm	2.0 mm	2.4 mm	2.7 mm	3.0 mm	3.5 mm Short	3.5 mm Standard	3.5 mm Broad	4.5 mm	4.5 mm
8 mm	✓									
10 mm	✓	✓								
12 mm		✓	✓							
15 mm			✓	✓						
18 mm				✓	✓					
21 mm					✓	✓	✓			
24 mm						✓	✓			
27 mm								✓		
30 mm								✓	✓	
33 mm										✓
	1-4 kg	3-7 kg	8-12 kg	13-17 kg	15-20 kg	18-25 kg	26-39 kg	40-55 kg	56+ kg	
										

All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultants' recommendations. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.

These are general guidelines, not definitive values, based on common blade radii observed in surgeon experience.

arthrexvetsystems.com

© 2025-04 Arthrex, Inc. All rights reserved. vLR1-000199-en-US_A

