



OrthoLine™ Bridge Plate System

Size Reference Chart








	Product Image	Item Number	Max Width (mm)	Length (mm)	Max Thickness (mm)	Bridge Length (mm)
Bridge Plate, Ti*, 1.6 mm Broad		VAR-3116BBP-11H	5.5	60	1.7	27
		VAR-3116BBP-12H	5.5	69	1.7	36
		VAR-3116BBP-13H	5.5	78	1.7	45
		VAR-3116BBP-14H	5.5	87	1.7	54
Bridge Plate, Ti*, 2.0 mm		VAR-3120BP-11H	6.1	74	2.3	33
		VAR-3120BP-12H	6.1	85	2.3	44
		VAR-3120BP-13H	6.1	96	2.3	55
		VAR-3120BP-14H	6.1	107	2.3	66

*Ti = titanium



OrthoLine™ Bridge Plate System







Size Reference Chart

	Product Image	Item Number	Max Width (mm)	Length (mm)	Max Thickness (mm)	Bridge Length (mm)
Bridge Plate, Ti*, 2.4 mm		VAR-3124BP-11H	7.4	92.4	2.8	42
		VAR-3124BP-12H	7.4	106.4	2.8	56
		VAR-3124BP-13H	7.4	120.3	2.8	70
Bridge Plate, Ti*, 3.0 mm		VAR-3130BP-11H	9.4	109	3.5	51
		VAR-3130BP-12H	9.4	126	3.5	68
		VAR-3130BP-13H	9.4	143	3.5	85
		VAR-3130BP-14H	9.4	160	3.5	102

*Ti = titanium

OrthoLine™ Bridge Plate System

Size Reference Chart

	Product Image	Item Number	Max Width (mm)	Length (mm)	Max Thickness (mm)	Bridge Length (mm)
Bridge Plate, SS†, 3.5 mm		VAR-3035BP-11H	12	128.6	4.0	60
		VAR-3035BP-12H	12	148.5	4.0	80
		VAR-3035BP-13H	12	168.6	4.0	100
		VAR-3035BP-14H	12	188.5	4.0	120
		VAR-3035BP-15H	12	208.6	4.0	140
		VAR-3035BP-16H	12	228.6	4.0	160

†SS = stainless steel