











# OrthoLine™ Fracture Plate System





## Size Chart for Acetabular Fracture Plates

	Product Image	Item Number	Length (mm)	Max Thickness (mm)	Max Width (mm)
Acetabular Fracture Plate, Ti, 1.6 mm		VAR- 3116AC-L	25	1.4	4.9
		VAR-3116AC-R			
Acetabular Fracture Plate, Ti, 2.0 mm		VAR-3120AC-L	31.9	1.5	6.2
		VAR-3120AC-R			
Acetabular Fracture Plate, Short, Ti, 2.4 mm		VAR- 3124SAC-L	33	1.9	7.6
		VAR-3124SAC-R			
Acetabular Fracture Plate, Ti, 2.4 mm		VAR-3124AC-L	40	1.9	7.6
		VAR-3124AC-R			



# OrthoLine™ Fracture Plate System

## Size Chart for Acetabular Fracture Plates

	Product Image	Item Number	Length (mm)	Max Thickness (mm)	Max Width (mm)
Acetabular Fracture Plate, Ti, 3.0 mm		VAR-3130AC-L	44	2.5	8.8
		VAR-3130AC-R			
Acetabular Fracture Plate, Broad, Ti, 3.0 mm		VAR-3130BAC-L	49.9	2.7	8.8
		VAR-3130BAC-R			

\*Ti = titanium