

OrthoLine™ Fracture Management System

Weight Reference Chart: Bridge Plates

	 1-4 kg	 4-8 kg	 8-12 kg	 13-19 kg	 20-30 kg	 31-45 kg	 45-55+ kg
Radius	1.4 mm, 1.6 mm	2.0 mm	2.0 mm, 2.4 mm	3.0 mm	3.0 mm, 3.5 mm	3.5 mm	3.5 mm
Humerus	1.4 mm, 1.6 mm	2.0 mm, 2.4 mm	2.4 mm, 3.0 mm	2.4 mm, 3.0 mm	3.5 mm	3.5 mm	
Femur	1.4 mm, 1.6 mm	2.0 mm, 2.4 mm	2.4 mm, 3.0 mm	2.4 mm, 3.0 mm	3.5 mm	3.5 mm	
Tibia	1.4 mm, 1.6 mm	2.0 mm, 2.4 mm	2.4 mm	2.4 mm	3.5 mm	3.5 mm	

All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultant's recommendation. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.