OrthoLine™ Pantarsal Arthrodesis Locking Plates with BacStop™ Antimicrobial Coating

Size Chart for Pantarsal Arthrodesis Plates

	Product Image	Item Number	Length (mm)	Maximum Head Thickness (mm)	Maximum Shaft Thickness (mm)	Maximum Head Width (mm)	Maximum Shaft Width (mm)	Screw sizes (mm)
Pantarsal Arthrodesis Plate, Ti, 1.6 mm		VAR-3116PT-L	60.8	3.7	3.7	9.1	6.1	1.6/1.4
		VAR-3116PT-R						
Pantarsal Arthrodesis Plate, Ti, 2.0 mm		VAR-3120PT-L	64.2	3.7	3.7	9.2	6.1	2.0/1.4
		VAR-3120PT-R						

Arthrex® Vet Systems

OrthoLine™ Pantarsal Arthrodesis Locking Plates with BacStop™ Antimicrobial Coating

Size Chart for Pantarsal Arthrodesis Plates


	Product Image	Item Number	Length (mm)	Maximum Head Thickness (mm)	Maximum Shaft Thickness (mm)	Maximum Head Width (mm)	Maximum Shaft Width (mm)	Screw sizes (mm)
Pantarsal Arthrodesis Plate, Ti, 2.4 mm		VAR-3124PT-L	72.6	4.0	4.0	10.5	7.4	2.4/2.0
		VAR-3124PT-R						
Pantarsal Arthrodesis Plate, Small, Ti, 3.0 mm		VAR-3130SMPT-L	85	4.8	4.8	13	9.4	3.0/2.0
		VAR-3130SMPT-R						
Pantarsal Arthrodesis Plate, Standard, Ti, 3.0 mm		VAR-3130SPT-L	88.3	4.8	4.8	14	9.4	3.0/2.4
		VAR-3130SPT-R						


OrthoLine™ Pantarsal Arthrodesis Locking Plates with BacStop™ Antimicrobial Coating

Size Chart for Pantarsal Arthrodesis Plates

	Product Image	Item Number	Length (mm)	Maximum Head Thickness (mm)	Maximum Shaft Thickness (mm)	Maximum Head Width (mm)	Maximum Shaft Width (mm)	Screw sizes (mm)
Pantarsal Arthrodesis Plate, Ti, 3.5 mm		VAR-3135PT-L	109.7	6.6	6.6	17.7	11.9	3.5/2.4
		VAR-3135PT-R						
Pantarsal Arthrodesis Plate, Broad, Ti, 3.5 mm		VAR-3135BPT-L	134.1	8.2	8.2	21.3	14.3	3.5/2.4
		VAR-3135BPT-R						

*Ti = titanium





© 2026-10 Arthrex, Inc. All rights reserved. vLR1-000219-en-US_A