
OrthoLine™ Fracture Management System

Weight Reference Chart: Acetabular Fracture Plate



1-4 kg	4-10 kg	8-18 kg	16-30 kg	25-39 kg	35-50 kg	45+ kg
1.4 mm	1.6 mm	2.0 mm	2.4 mm, short	2.4 mm	3.0 mm	3.0 mm, broad



All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultants' recommendations. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.

arthrexvetsystems.com

© 2026-07 Arthrex, Inc. All rights reserved. vLR1-000213-en-US_B

Arthrex® Vet Systems