
OrthoLine™ Arthrodesis Locking Plates With BacStop™ Antimicrobial Coating

Size Chart for Pancarpal Arthrodesis Plates

	Product Image	Item Number	Length (mm)	Maximum Head Thickness (mm)	Maximum Shaft Thickness (mm)	Maximum Head Width (mm)	Maximum Shaft Width (mm)
Pancarpal Arthrodesis Plate, Ti, 1.6 mm/1.4 mm		VAR-3116PC-L	49.6	1.7	2.1	8.3	6.0
		VAR-3116PC-R					
Pancarpal Arthrodesis Plate, Ti, 2.0 mm/1.4 mm		VAR-3120PC-L	55.6	1.9	2.2	9.1	6.0
		VAR-3120PC-R					
Pancarpal Arthrodesis Plate, Ti, 2.4 mm/2.0 mm		VAR-3124PC-L	78.8	2.6	2.6	11.2	7.4
		VAR-3124PC-R					
Pancarpal Arthrodesis Plate, Ti, 3.0 mm/2.0 mm		VAR-3130PC-L	88.8	2.9	3.9	14.1	9.4
		VAR-3130PC-R					



Arthrex® Vet Systems

OrthoLine™ Arthrodesis Locking Plates With BacStop™ Antimicrobial Coating

Size Chart for Pancarpal Arthrodesis Plates

	Product Image	Item Number	Length (mm)	Maximum Head Thickness (mm)	Maximum Shaft Thickness (mm)	Maximum Head Width (mm)	Maximum Shaft Width (mm)
Pancarpal Arthrodesis Plate, Ti, 3.5 mm/2.4 mm		VAR-3135PC-L	113.9	3.7	4.1	17.5	11.9
		VAR-3135PC-R					
Pancarpal Arthrodesis Plate, Broad, Ti, 3.5 mm/2.4 mm		VAR-3135BPC-L	137.8	4.5	5.0	21	14.3
		VAR-3135BPC-R					

*Ti = titanium

© 2026-10 Arthrex, Inc. All rights reserved. vLR1-000218-en-US_A