OrthoLine™ Arthrodesis Plating System With BacStop™ Antimicrobial Coating

Weight Reference Chart: Pancarpal Arthrodesis Plate

2-5 kg	3-10 kg	8-15 kg	15-25 kg	23-40 kg	40+ kg
1.6 mm/1.4 mm	2.0 mm/1.4 mm	2.4 mm/2.0 mm	3.5 mm/2.4 mm	3.5/2.4 mm	3.5 mm/2.4 mm, broad



All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultants' recommendations. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.

arthrexvetsystems.com

© 2026-10 Arthrex, Inc. All rights reserved. vLR1-000224-en-US_A

Arthrex® Vet Systems